

Store at
-20°C

#13859

Cellular Glutathione Detection Assay Kit

1 Kit (200 assays (96 well format))



Support: +1-978-867-2388 (U.S.)
cellsignal.com/support

Orders: 877-616-2355 (U.S.)
orders@cellsignal.com

For Research Use Only. Not for Use in Diagnostic Procedures.

Description: The Cellular Glutathione Detection Assay Kit employs the cell permeable dye monochlorobimane (MCB) to detect reduced glutathione (GSH) in cellular assays. MCB displays a high affinity for reduced glutathione and exhibits a very low fluorescent yield when free in solution. Upon binding to GSH, the dye exhibits a strong blue fluorescence that can be measured at an excitation wavelength of 380 nm and an emission wavelength of 460 nm. Fluorescent intensity correlates with sample GSH level. This kit can be used to either label cells directly or to detect GSH level in cell extracts. The assay can be easily applied in high-throughput plate, flow cytometry, or fluorescent imaging.

Specificity/Sensitivity: The Cellular Glutathione Detection Assay Kit is expected to detect the reduced glutathione in cellular assays across all species.

Background: The antioxidant glutathione is found in both reduced and oxidized states in cells. Reduced glutathione can play an important role in preventing cellular damage caused by reactive oxygen species, including free radicals and peroxides. Reduced glutathione (GSH) acts as an electron donor in the presence of free radicals and peroxides to become oxidized (GSSG). GSH also participates in redox signaling through the removal of the cellular second messenger H_2O_2 (1,2). Diminished glutathione levels are observed during the aging process and in oxidative stress-related diseases. The depletion of GSH is necessary for the progression of apoptosis that is mediated by various signaling pathways (3,4). Intracellular GSH levels can be a very useful indicator for overall cell health, proliferation, and death (2).


Background References:

- (1) Dickinson, D.A. and Forman, H.J. (2002) *Ann N Y Acad Sci* 973, 488-504.
- (2) Pompella, A. et al. (2003) *Biochem Pharmacol* 66, 1499-503.
- (3) Franco, R. et al. (2007) *J Biol Chem* 282, 30452-65.
- (4) Macho, A. et al. (1997) *J Immunol* 158, 4612-9.



Products Included	Item #	Quantity	Storage Temp
Reduced Glutathione Standard	13870	1 Ea	-20°C
Glutathione-S-Transferase	13882	1 Ea	-20°C
Tris Assay Buffer	13865	25 ml	4°C
Digitonin Lysis Buffer	13899	11 ml	-20°C
Monochlorobimane	13887	1 Ea	-20°C

Note: Upon receipt, Tris Assay Buffer (#13865) should be removed from #13859 and stored at 4°C. And remaining components should be stored at -20°C.

Note: All components in this kit are stable for 36 months when stored at the recommended temperature and left unused.



◀ Figure 1. Reduced Glutathione Standard was diluted in Tris Assay Buffer and samples were assayed to create a standard curve. This standard curve is for demonstration purposes only; in order to accurately determine glutathione levels users should generate a standard curve for each sample set.



◀ Figure 2. Raw 264.7, HeLa, and Jurkat cells were seeded in 96-well plates at various cell densities. Cellular glutathione levels were determined using the Cellular Glutathione Detection Assay Kit using both cell extracts (upper) and live cells (lower).

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

cellsignal.com




Figure 3. The relationship between the protein concentration of lysates from HeLa Cells, untreated and treated with Staurosporine #9953 (2.0 μ M), and RFU as determined by the Cellular Glutathione Detection Assay Kit is shown.




Figure 4. HeLa cells (4×10^4 cells/well) were treated with terfenadine (indicated concentrations, 4 hr) and reduced glutathione levels were determined using the Cellular Glutathione Detection Assay Kit. Relative fluorescence units (RFU) were measured using a fluorescent plate reader with an excitation wavelength of 380 nm and an emission wavelength of about 460 nm.




Figure 5. Jurkat cells, untreated (green) or treated with terfenadine (200 μ M, 4 hr; blue), were labeled with MCB (100 μ M, 30 min) and analyzed using direct flow cytometry at an excitation wavelength of 380 nm and an emission wavelength of about 460 nm.




Figure 6. Jurkat cells, untreated or terfenadine-treated (200 μ M, 4 hr), were labeled with MCB (100 μ M, 30 min) and Propidium Iodide (PI)/RNase Staining Solution #4087 and analyzed under a fluorescent microscope. Live cells (untreated) show high reduced-glutathione level (blue) while dead cells (treated) stained with PI (red) due to their compromised plasma

Cellular Glutathione Detection Assay Kit Protocol

A Solutions and Reagents

- Reduced Glutathione Standard:** Reconstitute lyophilized product in 100 μ l of Tris Assay Buffer to create a 100 mM solution.
- Glutathione-S-Transferase:** Reconstitute lyophilized product with 110 μ l of Tris Assay Buffer.
- Tris Assay Buffer:** Bring to room temperature before use.
- Digitonin Lysis Buffer:** Thaw before use, and keep on ice during use.
- Monochlorobimane:** Reconstitute lyophilized product with 44 μ l of DMSO
- Create a **Reduced Glutathione Working Solution** as follows (see Table 1).

TABLE 1: GSH Working Solution for one 96-well plate (100 assays)

Kit Component	Volume (μ l)
Monochlorobimane	20
Glutathione-S-Transferase	50
Tris Assay Buffer	4930
Total Volume (μl)	5000

NOTE: To create a GSH standard curve, prepare a 100 μ M GSH Standard Solution by diluting the 100 mM Reduced Glutathione Standard 1:1000 in Tris Assay Buffer. Run the standard curve in triplicate using 50 μ l of standard per well. Create a five point standard curve by starting at 100 μ M and using a 1:2 serial dilution in Tris Assay Buffer. Include a blank control. The GSH Standard Solution is not stable at this concentration and should be used within 2 hours or discarded.

Additional Reagents (not supplied): DMSO (#12611)

B Preparing Cell Extracts

- Preparation of cell extracts (10 cm dish)**
 - Apply desired treatment protocol and prepare a control sample without treatment.
 - Following treatment, aspirate media and add 1.0 ml Digitonin Lysis Buffer to each sample, adjusting Digitonin Lysis Buffer volume when using a small surface area vessel.

NOTE: If using suspended or detached cells, collect and centrifuge the cell suspension for 5 min at 1,200 rpm. Discard the supernatant and add Digitonin Lysis Buffer to the cell pellet.
 - Incubate the cell lysates on ice for 5-10 min.
 - Scrape cells off the plate and transfer suspension to a microcentrifuge tube.
 - Centrifuge cell lysates for 10 min at 14,000 rpm at 4°C.
 - Collect the supernatant and run assays or store at -80°C for later use.

NOTE: Use a 10 kDa MWCO centrifugal filter and follow the manufacturer's instructions for cells that require an additional filtration step.
- Preparation of cell extracts from 96-well assay plate**
 - Following treatment, centrifuge assay plate for 10 min at 1,200 rpm and carefully remove the supernatant.
 - Add 50 μ l Digitonin Lysis Buffer to each well and incubate at room temperature on shaker for 15 min.
 - Transfer samples into microcentrifuge tubes and centrifuge samples for 10 min at 14,000 rpm at 4°C.
 - Collect the supernatants and immediately proceed to GSH assay.

C Test Procedure

- Detection of Reduced Glutathione in Cell Extracts**
 - Mix Tris Assay Buffer with test sample, GSH standard, and Working Solution in a black 96-well plate with a clear bottom to final volume of 100 μ l (See Table 2).

TABLE 2: Detection and Standard Curve Assays

Sample	Blank	Standard	Sample
Test Sample (μl)	—	—	1-50
GSH Standard (μl)	—	50	—
Working Solution (μl)	50	50	50
Assay Buffer (μl)	50	—	50 minus sample volume
Total Volume (μl)	100	100	100

- Incubate the assay plate at room temperature for 60 min, protected from light.
- Read plate with a plate reader at an excitation wavelength around 380 nm and emission wavelength of 485 nm.

2. Detection of Reduced Glutathione (GSH) in Live Cells

The following procedure measures reduced glutathione in live cells using a 96-well plate. The same protocol can be modified for fluorescent imaging and flow cytometry detection.

- Culture adherent or suspension cells in a 96-well black plate with clear bottom in duplicate for minimum of 16 hr. Determine the background reading by using at least 2 wells containing growth media without cells.

NOTE: For best results, a cell number titration is recommended to determine the optimal cell seeding density. We recommend not exceeding 1×10^5 cells per well.
- Create a Live Cell Staining Working Solution by diluting the reconstituted Monochlorobimane stock solution 1:50 in Tris Assay Buffer. For example, dilute 20 μ l of Monochlorobimane stock solution to 1 ml Tris Assay Buffer to make 1 ml Live Cell Staining Working Solution, which is enough for one 96-well plate. Mix reagents by gently tapping each side of the plate.
- Add 10 μ l of Live Cell Staining Working Solution to each well containing 100 μ l of growth media.
- Incubate samples in a CO₂ incubator for a minimum of 30 min. For best results, measure the fluorescence intensity every 30 min for up to 3 hr or until the signal has reached a plateau.