

Serine/Threonine Kinase Substrate Screening Kit

☑ 2 assays

(After initial use, store plates at appropriate temperature)



Cell Signaling
TECHNOLOGY®

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 10/04/07

This product is for *in vitro* research use only and is not intended for use in humans or animals.

| Products Included | Products # | Kit Quantity |
|---|------------|---------------------|
| Antibody Plates: 96-well plates | | 6 |
| 87 Phospho-Specific Antibodies | See table | 2 X 100 µl each |
| Peptide plates: 96-well plates | 9804 | 3 |
| 87 Phospho-Peptides (Biotinylated) | | 50 µl, 12.5 µM each |
| 87 X 2 Nonphospho-Peptides (Biotinylated) | | 50 µl, 12.5 µM each |
| Kinase Buffer (10X) | 9802 | 15 ml |
| ATP (10 mM) | 9804 | 1 ml |
| CD Rom; data, overview | | 1 |
| Kinase of interest (not included) | | |

Background: In order to expedite the isolation of substrates for novel kinases, Cell Signaling Technology has developed a kinase/substrate development kit. Numerous serine/threonine kinases have become pharmaceutically important, including GSK, Akt and p38. With the cloning of hundreds of uncharacterized serine/threonine kinases, drug discovery researchers are seeking to quickly find substrates for high throughput screening programs.

With a total of 87 different peptide/phospho-specific antibody pairs on the array, a very diverse group of substrates is represented. Signal-to-noise ratios have been optimized and range from 10 to over 1000. The phospho-specific antibody and peptide substrate pair can be used for subsequent high throughput peptide based kinase screens (e.g., FP, TRF, AlphaScreen).


Description: The Serine/Threonine Kinase Substrate Screening Kit provides 87 biotinylated nonphospho-peptides as kinase substrates. Also included are biotinylated phospho and nonphospho-peptides of identical sequence as positive and negative controls respectively. Potential substrate phosphorylation is detected by the corresponding phospho-specific antibodies.

The paired phospho-specific antibody and peptide substrate can be used for high throughput peptide based kinase screens (e.g. FP, TRF, AlphaScreen).

Specificity/Sensitivity: 10 to 1000+ fold signal-to-noise ratio of phospho- versus nonphospho-peptides measured by DELFIA® assay.

Source/Purification: Polyclonal antibodies are produced by immunizing rabbits with synthetic phospho-peptides (KLH-coupled) and are purified by protein A and peptide affinity chromatography. Monoclonal antibodies are derived from mice immunized with synthetic phospho-peptides coupled to KLH.

Quality Control: The quality of the peptides was evaluated by reverse-phase HPLC and by mass spectrometry.



Kinase of interest is incubated with biotinylated substrate peptides. Phosphorylated products are captured on streptavidin-coated plates. Cell Signaling Technology antibodies serve as primary detection reagents; DELFIA® Europium (Eu) anti-rabbit or anti-mouse chelates are used as secondary antibodies. Positive hits from Serine/Threonine Kinase Substrate Screening Kit (#7400) can be used in any subsequent peptide based kinase assay.

Storage: Antibodies are supplied in 50 mM Tris-HCl (pH 7.4), 150 mM NaCl (TBS), 0.05% Tween-20, 1.0 µg/ml BSA. Store at 4°C. Do not aliquot the antibodies.

Peptides are supplied at 12.5 µM in 1X kinase buffer (25 mM Tris (pH 7.5), 5 mM beta-glycerolphosphate, 2 mM DTT, 0.1 mM Na₃VO₄, 5 mM MgCl₂). Store peptides at -20°C.

**KEY:**

Plate, Row, and Column designation refer to positions on Kinase Substrate Screening Kit

P = CST phospho-peptide number (positive control found in wells 1, 4, 7, 10)**NP** = CST nonphospho-peptide number (negative control wells found in wells 2, 5, 8, 11 and kinase substrates found in wells 3, 6, 9, 12)**SNR** = Signal-to-noise ratio (P DELFIA® units/NP DELFIA® units)

| Plate | Row | Column | P# | NP# | CST Ab # | Target/Antibody Description | S/N |
|-------|-----|----------|------|------|----------|--|------|
| 1 | a | 1,2,3 | 872 | 870 | 2974 | Phospho-mTOR (Ser2481) Antibody | 1343 |
| 1 | a | 4,5,6 | 910 | 907 | 2301 | Phospho-DARPP-32 (Thr75) Antibody | 394 |
| 1 | a | 7,8,9 | 826 | 821 | 3081 | Phospho-C/EBP β (Ser105) Antibody (Rat Specific) | 160 |
| 1 | a | 10,11,12 | 848 | 844 | 9291 | Phospho-Bad (Ser112) Antibody | 78 |
| 1 | b | 1,2,3 | 1132 | 1131 | 3361 | Phospho-CaMKII (Thr286) Antibody | 948 |
| 1 | b | 4,5,6 | 911 | 908 | 3141 | Phospho-Ezrin (Thr567)/Radixin (Thr564)/Moesin (Thr558) Antibody | 340 |
| 1 | b | 7,8,9 | 1166 | 1133 | 9336 | Phospho-GSK-3 β (Ser9) Antibody | 143 |
| 1 | b | 10,11,12 | 850 | 846 | 2341 | Phospho-Chk1 (Ser345) Antibody | 76 |
| 1 | c | 1,2,3 | 842 | 836 | 9121 | Phospho-MEK1/2 (Ser217/221) Antibody | 840 |
| 1 | c | 4,5,6 | 871 | 869 | 9571 | Phospho-eNOS (Ser1177) Antibody | 290 |
| 1 | c | 7,8,9 | 489 | 485 | 9601 | Phospho-(Ser) 14-3-3 Binding Motif Antibody | 119 |
| 1 | c | 10,11,12 | 1152 | 1151 | 3114 | Phospho-VASP (Ser239) Antibody | 74 |
| 1 | d | 1,2,3 | 1168 | 1167 | 9376 | Phospho-PKC δ/θ (Ser643/676) Antibody | 698 |
| 1 | d | 4,5,6 | 489 | 485 | 9528 | Phospho-cdc25C (Ser216) Antibody | 243 |
| 1 | d | 7,8,9 | 503 | 505 | 9611 | Phospho-(Ser/Thr) Akt Substrate Antibody | 117 |
| 1 | d | 10,11,12 | 829 | 824 | 9461 | Phospho-FoxO1 (Ser256) Antibody | 72 |
| 1 | e | 1,2,3 | 1148 | 1147 | 1932 | Phospho-PLK (Ser137) Antibody | 561 |
| 1 | e | 4,5,6 | 888 | 890 | 9631 | Phospho-(Ser/Thr) Phe Antibody | 228 |
| 1 | e | 7,8,9 | 1128 | 1127 | 9591 | Phospho-MSK1 (Ser376) Antibody | 97 |
| 1 | e | 10,11,12 | 1150 | 1149 | 9151 | Phospho-SEK1/MKK4 (Thr261) Antibody | 71 |
| 1 | f | 1,2,3 | 861 | 865 | 9297 | Phospho-Bad (Ser155) Antibody | 533 |
| 1 | f | 4,5,6 | | | | | |
| 1 | f | 7,8,9 | 922 | 924 | 9287 | Phospho-p53 (Ser20) Antibody | 95 |
| 1 | f | 10,11,12 | 493 | 494 | 9131 | Phospho-Stat3 (Tyr705) Antibody | 69 |
| 1 | g | 1,2,3 | 945 | 940 | 9455 | Phospho-4E-BP1 (Thr70) Antibody | 512 |
| 1 | g | 4,5,6 | 1126 | 1125 | 9181 | Phospho-Elk-1 (Ser383) Antibody | 194 |
| 1 | g | 7,8,9 | 828 | 823 | 2441 | Phospho-eIF4G (Ser1108) Antibody | 91 |
| 1 | g | 10,11,12 | 851 | 847 | 9231 | Phospho-MKK3/MKK6 (Ser189/207) Antibody | 65 |
| 1 | h | 1,2,3 | 1130 | 1129 | 2864 | Phospho-c-Abl (Thr735) Antibody | 463 |
| 1 | h | 4,5,6 | 923 | 915 | 9307 | Phospho-Rb (Ser780) Antibody | 170 |
| 1 | h | 7,8,9 | 946 | 943 | 9114 | Phospho-cdc2 (Thr161) Antibody | 84 |
| 1 | h | 10,11,12 | 862 | 858 | 9331 | Phospho-GSK-3 α/β (Ser21/9) Antibody | 61 |
| 2 | a | 1,2,3 | 1215 | 1136 | 9204 | Phospho-p70 S6 Kinase (Thr421/Ser424) Antibody | 56 |
| 2 | a | 4,5,6 | 1128 | 1127 | 9631 | Phospho-(Ser/Thr) Phe Antibody | 43 |
| 2 | a | 7,8,9 | 909 | 905 | 2531 | Phospho-AMPK- α (Thr172) Antibody | 28 |
| 2 | a | 10,11,12 | 488 | 484 | 2291 | Phospho-(Thr) PDK1 Substrate Antibody | 16 |

KEY:

Plate, Row, and Column designation refer to positions on Kinase Substrate Screening Kit

P = CST phospho-peptide number (positive control found in wells 1, 4, 7, 10)

NP = CST nonphospho-peptide number (negative control wells found in wells 2, 5, 8, 11 and kinase substrates found in wells 3, 6, 9, 12)

SNR = Signal-to-noise ratio (P DELFIA® units/NP DELFIA® units)


| Plate | Row | Column | P# | NP# | CST Ab # | Target/Antibody Description | S/N |
|-------|-----|----------|------|------|----------|--|-----|
| 2 | b | 1,2,3 | 1144 | 1143 | 3068 | Phospho-Aurora A (Thr288)/Aurora B (Thr232)/Aurora C (Thr198) Antibody | 156 |
| 2 | b | 4,5,6 | 1148 | 1147 | 2261 | Phospho-(Ser) PKC Substrate Antibody | 40 |
| 2 | b | 7,8,9 | 825 | 820 | 9295 | Phospho-Bad (Ser136) Antibody | 26 |
| 2 | b | 10,11,12 | 863 | 859 | 2401 | Phospho-HSP27 (Ser82) Antibody | 58 |
| 2 | c | 1,2,3 | 503 | 505 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 20 |
| 2 | c | 4,5,6 | 849 | 845 | 9164 | Phospho-c-Jun (Ser73) Antibody | 36 |
| 2 | c | 7,8,9 | 1182 | 1181 | 9101 | Phospho-p44/42 MAP Kinase (Thr202/Tyr204) Antibody | 25 |
| 2 | c | 10,11,12 | 1396 | 1395 | 2981 | Phospho-(Ser) Arg-X-Tyr/Phe-X-pSer Motif Antibody | 13 |
| 2 | d | 1,2,3 | 490 | 498 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 53 |
| 2 | d | 4,5,6 | 827 | 822 | 3311 | Phospho-Cofilin (Ser3) Antibody | 34 |
| 2 | d | 7,8,9 | 873 | 877 | 9551 | Phospho-PTEN (Ser380) Antibody | 22 |
| 2 | d | 10,11,12 | 1178 | 1177 | 9464 | Phospho-FoxO1 (Thr24)/FoxO3a (Thr32) Antibody | 9 |
| 2 | e | 1,2,3 | 888 | 890 | 9341 | Phospho-p90RSK (Ser380) Antibody | 50 |
| 2 | e | 4,5,6 | 1146 | 1145 | 2661 | Phospho-Chk2 (Thr68) Antibody | 30 |
| 2 | e | 7,8,9 | 1180 | 1179 | 9211 | Phospho-p38 MAP Kinase (Thr180/Tyr182) Antibody | 22 |
| 2 | e | 10,11,12 | 1032 | 1031 | 9221 | Phospho-ATF-2 (Thr71) Antibody | 6 |
| 2 | f | 1,2,3 | 491 | 500 | 2261 | Phospho-(Ser) PKC Substrate Antibody | 49 |
| 2 | f | 4,5,6 | 1124 | 1123 | 9261 | Phospho-c-Jun (Ser63) II Antibody | 28 |
| 2 | f | 7,8,9 | 1138 | 1137 | 9346 | Phospho-p90RSK (Thr573) Antibody | 10 |
| 2 | f | 10,11,12 | 884 | 880 | 9741 | Phospho-eIF4E (Ser209) Antibody | 33 |
| 2 | g | 1,2,3 | 1140 | 1139 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 45 |
| 2 | g | 4,5,6 | 1184 | 1183 | 9251 | Phospho-SAPK/JNK (Thr183/Tyr185) Antibody | 28 |
| 2 | g | 7,8,9 | 1132 | 1131 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 267 |
| 2 | g | 10,11,12 | 491 | 500 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 20 |
| 2 | h | 1,2,3 | 1214 | 1213 | 2851 | Phospho-(Ser/Thr) ATM/ATR Substrate Antibody | 203 |
| 2 | h | 4,5,6 | 853 | 855 | 9631 | Phospho-(Ser/Thr) Phe Antibody | 20 |
| 2 | h | 7,8,9 | 864 | 876 | 2851 | Phospho-(Ser/Thr) ATM/ATR Substrate Antibody | 24 |
| 2 | h | 10,11,12 | 874 | 878 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 32 |
| 3 | a | 1,2,3 | 875 | 879 | 9621 | Phospho-(Ser/Thr) PKA Substrate Antibody | 76 |
| 3 | a | 4,5,6 | 1170 | 1169 | 2851 | Phospho-(Ser/Thr) ATM/ATR Substrate Antibody | 140 |
| 3 | a | 7,8,9 | 1176 | 1175 | 9631 | Phospho-(Ser/Thr) Phe Antibody | 20 |
| 3 | a | 10,11,12 | | | | intentionally left blank | |
| 3 | b | 1,2,3 | | | | intentionally left blank | |
| 3 | b | 4,5,6 | | | | intentionally left blank | |
| 3 | b | 7,8,9 | | | | intentionally left blank | |
| 3 | b | 10,11,12 | | | | intentionally left blank | |
| 3 | c | 1,2,3 | | | | intentionally left blank | |

**KEY:**

Plate, Row, and Column designation refer to positions on Kinase Substrate Screening Kit

P = CST phospho-peptide number (positive control found in wells 1, 4, 7, 10)**NP** = CST nonphospho-peptide number (negative control wells found in wells 2, 5, 8, 11 and kinase substrates found in wells 3, 6, 9, 12)**SNR** = Signal-to-noise ratio (P DELFIA® units/NP DELFIA® units)

| Plate | Row | Column | P# | NP# | CST Ab # | Target/Antibody Description | S/N |
|-------|-----|----------|------|------|----------|---|-----|
| 3 | c | 4,5,6 | | | | intentionally left blank | |
| 3 | c | 7,8,9 | | | | intentionally left blank | |
| 3 | c | 10,11,12 | | | | intentionally left blank | |
| 3 | d | 1,2,3 | 864 | 876 | 9286 | Phospho-p53 (Ser15) (16G8) mAb | 107 |
| 3 | d | 4,5,6 | 1184 | 1183 | 9255 | Phospho-SAPK/JNK (Thr183/Tyr185) (G9) mAb | 59 |
| 3 | d | 7,8,9 | 1130 | 1129 | 9386 | Phospho-Threonine (42H4) mAb | 25 |
| 3 | d | 10,11,12 | 487 | 483 | 9391 | Phospho-Threonine-Proline mAb (P-Thr-Pro-101) | 47 |
| 3 | e | 1,2,3 | 839 | 833 | 9246 | Phospho-I κ B- α (Ser32/36) (5A5) mAb | 83 |
| 3 | e | 4,5,6 | 1180 | 1179 | 9216 | Phospho-p38 MAPK (Thr180/Tyr182) (28B10) mAb | 36 |
| 3 | e | 7,8,9 | 1135 | 1134 | 9386 | Phospho-Threonine (42H4) mAb | 20 |
| 3 | e | 10,11,12 | 825 | 820 | 9606 | Phospho-(Ser) 14-3-3 Binding Motif (4E2) mAb | 11 |
| 3 | f | 1,2,3 | 838 | 832 | 9706 | Phospho-Histone H3 (Ser10) (6G3) mAb | 83 |
| 3 | f | 4,5,6 | 912 | 913 | 9634 | Phospho-(Ser/Thr) PDK1 Docking Motif (18A2) mAb | 17 |
| 3 | f | 7,8,9 | 1032 | 1031 | 9391 | Phospho-Threonine-Proline mAb (P-Thr-Pro-101) | 17 |
| 3 | f | 10,11,12 | 1132 | 1131 | 9391 | Phospho-Threonine-Proline mAb (P-Thr-Pro-101) | 24 |
| 3 | g | 1,2,3 | 489 | 485 | 9606 | Phospho-(Ser) 14-3-3 Binding Motif (4E2) mAb | 62 |
| 3 | g | 4,5,6 | 488 | 484 | 9386 | Phospho-Threonine (42H4) mAb | 31 |
| 3 | g | 7,8,9 | 1184 | 1183 | 9391 | Phospho-Threonine-Proline mAb (P-Thr-Pro-101) | 16 |
| 3 | g | 10,11,12 | 888 | 890 | 9634 | Phospho-(Ser/Thr) PDK1 Docking Motif (18A2) mAb | 8 |
| 3 | h | 1,2,3 | 1182 | 1181 | 9106 | Phospho-p44/42 MAPK (Thr202/Tyr204) (E10) mAb | 60 |
| 3 | h | 4,5,6 | 1146 | 1145 | 9386 | Phospho-Threonine (42H4) mAb | 31 |
| 3 | h | 7,8,9 | 1174 | 1173 | 9136 | Phospho-Stat3 (Ser727) (6E4) mAb | 15 |
| 3 | h | 10,11,12 | 1135 | 1134 | 9391 | Phospho-Threonine-Proline Mouse mAb (P-Thr-Pro-101) | 25 |



Protocol for Serine/Threonine Kinase Substrate Screening Kit

Suggested Kinase Assay Conditions:

Combine 10 μ M peptide substrate, 200 μ M ATP and 10-100 nM Kinase in a 25 μ l reaction. Incubate at 30°C for 30 minutes.

Peptide substrate is supplied in 1X kinase buffer, which was diluted from 10X Kinase Buffer #9802.

Note: Optimal incubation times and enzyme concentrations must be determined empirically for each particular kinase.

A Additional Solutions and Reagents (Not included)

- 1. Wash Buffer:** 1X TBS, 0.05% Tween-20 (TBS/T)
- Bovine Serum Albumin (BSA)
- 3. Stop Buffer:** tBST, 10% BSA, 100 mM EDTA
- DELFLIA® Europium-labeled Anti-rabbit IgG (PerkinElmer Life Sciences #AD0105)
- DELFLIA® Europium-labeled Anti-mouse IgG (PerkinElmer Life Sciences #AD0124)
- DELFLIA® Enhancement Solution (PerkinElmer Life Sciences #1244-105)
- DELFLIA® Streptavidin coated, 96-well, yellow plate (PerkinElmer Life Sciences AAAND-0005)

DELFLIA® is a registered trademark of PerkinElmer Life Sciences

B Protocol

1. Set up Reaction cocktails (Mock/No Kinase and Plus Kinase) in tubes on ice.

For each 25 μ l reaction, add the following to the Reaction cocktails:

0.5 μ l ATP (10 mM)

10-100 nM Kinase (Plus Kinase Cocktail Only)

Each 96-well plate will require Plus Kinase Reaction cocktail for 32 x 5 μ l reactions and Mock/No Kinase Reaction cocktail for 64 x 5 μ l. It is recommended that sufficient Reaction cocktail solution be prepared for 33 and 65 reactions respectively.

2. *Prepare substrate plates on ice. Spin peptide plates 5 minutes at 1000 rpm.

For kinase reaction with 10 μ M peptide substrate, transfer 20 μ l peptide solution (12.5 μ M) to two 96-well plates on ice. Maintain primary plate configuration with all transfers.

3. Add 5 μ l Reaction cocktails to peptide substrate plates. Keep on ice. Add Mock/No Kinase reaction cocktail to plate columns 1-2, 4-5, 7-8, 10-11. Add Plus Kinase Reaction cocktail to plate columns 3, 6, 9, 12. Mix with gentle agitation or pipetting.
4. Incubate reaction plate at 30°C for 30 minutes.
5. *Spin plates 5 minutes at 1000 rpm.
6. Add 100 μ l Stop Buffer to each well and mix by pipetting.
7. Transfer 5-100 μ l of each reaction to 96-well streptavidin-coated plate and incubate at room temperature for 30 minutes.

Indicated DELFLIA® signal-to-noise values were obtained by capturing 5 pmoles peptide on streptavidin. DELFLIA® signal-to-noise ratio was calculated as the ratio between phospho-peptide and nonphospho-peptide DELFLIA® signals.

8. *Wash three times with 100 μ l/well TBS/T.
9. Add 100 μ l primary antibody from antibody plate to matched well of streptavidin-coated plate.
10. Mix and incubate at room temperature with rocking for 60 minutes.
11. *Wash three times with 100 μ l/well TBS/T.
12. Add 100 μ l/well diluted Europium-labeled anti-rabbit or anti-mouse IgG antibody. Dilute the anti-rabbit IgG antibody 1:1000 and the anti-mouse IgG antibody 1:500 in DELFLIA® Assay buffer.
13. Incubate at room temperature for 30 minutes.
14. *Wash five times with 100 μ l/well TBS/T.
15. Add 100 μ l/well DELFLIA® Enhancement Solution.
16. Incubate at room temperature for 5 minutes.
17. Detect 615 nm fluorescence emission with appropriate Time-Resolved Plate Reader.

***IMPORTANT: Use of an automated microplate washer as well as centrifugation of plates when appropriate, greatly improves reproducibility.**

Please contact Cell Signaling Technology for HTS-ready antibodies (PBS formulated and carrier-free), and detailed peptide substrate sequence information.
Email: drugdiscovery@cellsignal.com