

Store at  
-20°C

# SignalScan™ Peptide Mix (SARS-CoV-2)



#23733

96 µl  
(96 assays)Support: +1-978-867-2388 (U.S.)  
www.cellsignal.com/supportOrders: 877-616-2355 (U.S.)  
orders@cellsignal.com

rev 04/09/21

## For Research Use Only. Not For Use In Diagnostic Procedures.

**Description:** The SignalScan™ Peptide Mix (SARS-CoV-2) enables targeted analysis of selected peptides in human-sourced cells or tissues. Peptides and phosphopeptides derived from SARS-CoV-2 spike and nucleocapsid proteins as well as nine human host response proteins are included. This mixture contains 96 pmol each of 25 heavy-labeled (<sup>13</sup>C, <sup>15</sup>N) Stable Isotope Standard (SIS) peptides, quantified by amino acid analysis (AAA). Inclusion of SIS peptides enhances the sensitivity of the LCMS assay by acting as a trigger for targeted MS/MS analysis. Targeted assays ensure that peptides present in the biological sample of interest are consistently observed and sensitivity is maximized. SIS-assisted scanning provides the sensitivity of a parallel reaction monitoring (PRM) assay without use of rigid retention time windows. Full peptide list, detailed protocol, and files for method set-up and data analysis are available for download at <http://media.cellsignal.com/www/zip/proteomics/signalscan-peptide-mix.zip>.

**Background:** The cause of the COVID-19 pandemic is a novel and highly pathogenic coronavirus, termed SARS-CoV-2 (severe acute respiratory syndrome coronavirus-2). SARS-CoV-2 is a member of the Coronaviridae family of viruses (1). The genome of SARS-CoV-2 is similar to other coronaviruses, and is comprised of four key structural proteins: S, the spike protein, E, the envelope protein, M, the membrane protein, and N, the nucleocapsid protein (2). Coronavirus spike proteins are class I fusion proteins and harbor an ectodomain, a trans-membrane domain, and an intracellular tail (3,4). The highly glycosylated ectodomain projects from the viral envelope surface and facilitates attachment and fusion with the host cell plasma membrane. The ectodomain can be further subdivided into host receptor-binding domain (RBD) (S1) and membrane-fusion (S2) subunits, which are produced upon proteolysis by host proteases at S1/S2 and S2' sites. S1 and S2 subunits remain associated after cleavage and assemble into crown-like homotrimers (2,4). In humans, both SARS-CoV and SARS-CoV-2 spike proteins utilize the angiotensin-converting enzyme 2 (ACE2) protein as a receptor for cellular entry (5-7). The cell's subsequent host response to SARS-CoV-2 infection involves complex signaling cascades, observable through perturbations in protein phosphorylation patterns (8,9) that lead to interferon and cytokine production (10).


Protein	PTM	Peptide
Spike		TQLPPAYTNSFTR
Spike		VGGNYNYLYR
Spike	R815 uncleaved	RSFIEDLLFNK
Spike		SFIEDLLFNK
Spike		LIANQFNSAIGK
Nucleocapsid		RPQGLPNNTASWFTALTQHGK
Nucleocapsid		ADETQALPQR
Nucleocapsid		GFYAEGSR
Nucleocapsid	pS176	GFYAEGS*R
Nucleocapsid	pT198	NST*PGSSR
Nucleocapsid	pS206	GTS*PAR
ACE2		GDYEVNGVDGYDYSR
ACE2		NSFVGWSTDWSPYADQSIK
IL6		EALAENNLNLPK
IL6		YILDGISALR
IL8		ELCLDPK
Interferon alpha		SLSLSTNLQER
Interferon beta		LLWQLNGR
Jak1	pY1034	AIETDKEY*YTVK
Jak1	pY1034/pY1035	AIETDKEY*Y*TVK
STAT1	pS727	LQTTDNLLPMS*PEEFDEVS
STAT1	pY701	GTGY*IKTELISVSEVHPSR
IRAK4	pT345	FAQTVMT*SR
IRAK4	pT345/pS346	FAQTVMT*S*R
PERK	pS715	EHIEIAPS*PQR

**Storage:** This product is stable for 12 months when stored at -20°C. Aliquot to avoid multiple freeze/thaw cycles.

### Background References:

- (1) Zhou, P. et al. (2020) *Nature* 579, 270-273.
- (2) Tortorici, M.A. and Velesler, D. (2019) *Adv Virus Res* 105, 93-116.
- (3) Li, F. et al. (2006) *J Virol* 80, 6794-800.
- (4) Li, F. (2016) *Annu Rev Virol* 3, 237-261.
- (5) Shang, J. et al. (2020) *Nature* 581, 221-224.
- (6) Wrapp, D. et al. (2020) *Science* 367, 1260-1263.
- (7) Yan, R. et al. (2020) *Science* 367, 1444-1448.
- (8) Davidson, A.D. et al. (2020) *Genome Med* 12, 68.
- (9) Bouhaddou, M. et al. (2020) *Cell* 182, 685-712. e19.
- (10) Blanco-Melo, D. et al. (2020) *Cell* 181, 1036-1045.e9.

◀ Peptides included in SignalScan™ Peptide Mix (SARS-CoV-2). Peptides in blue are derived from SARS-CoV-2 proteins. Peptides in black represent human proteins. All peptides are stable-isotope labeled on their carboxy-terminal amino acid (<sup>15</sup>N, <sup>13</sup>C<sub>6</sub>-Arg or <sup>15</sup>N, <sup>13</sup>C<sub>6</sub>-Lys), designated by bold **R** or **K**. Phosphorylated residues are designated with an asterisk (\*).



◀ SignalScan™ Workflow Steps. The SignalScan™ Peptide Mix (SARS-CoV-2) is compatible with any human sourced biological material. Extract proteins from the tissue and digest them into peptides using any standard trypsin-based protocol. Spike in 1 pmol of the heavy Stable Isotope Standard (SIS) peptides provided prior to LC-MS/MS analysis.

SureQuant is a registered trademark of Thermo Fisher Scientific Inc.

Thank you for your recent purchase. If you would like to provide a review visit [cellsignal.com/comments](http://cellsignal.com/comments).

[www.cellsignal.com](http://www.cellsignal.com)

© 2020 Cell Signaling Technology, Inc.

SignalScan and Cell Signaling Technology are trademarks of Cell Signaling Technology, Inc.

Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide Species Cross-Reactivity: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.

## SignalScan™ Peptide Mix (SARS-CoV-2) Protocol

### A Product Provided

SignalScan™ Peptide Mix (SARS-CoV-2) at a concentration of 1 pmol/μl per peptide in 0.1% trifluoroacetic acid, 25% acetonitrile. Peptides are labeled with <sup>13</sup>C, <sup>15</sup>N at carboxy-terminal lysine or arginine residues and are quantified by amino acid analysis.

### B Downloadable Analysis Package Contains

- Peptide list
- SignalScan™ Stable Isotope Standard Assisted Scanning (SISAS) protocol
- SureQuant method screenshots
- Skyline templates
- Survey Run\_Targeted Mass.csv
- SureQuant\_Targeted Mass\_Intensity.csv example files (four files: K2+, K3+, R2+, R3+)
- SureQuant\_Targeted Mass\_Trigger.csv (four files: K2+, K3+, R2+, R3+)

### C Sample Preparation and Analysis

- Dilute SignalScan™ Peptide Mix (SARS-CoV-2) at least 5x in your LC-MS loading buffer (not provided, 0.1% TFA recommended) prior to LC-MS injection to reduce the acetonitrile concentration.
- For SureQuant analysis, dilute 1 pmol SignalScan™ Peptide Mix (SARS-CoV-2) into 1-2 ug purified peptide sample (from tissue, cell lysate, or other sample digested with trypsin) for each LC-MS injection.
- For details of SureQuant method LCMS setup and data analysis, download analysis package from <http://media.cellsignal.com/www/zip/proteomics/signalscan-peptide-mix.zip> and follow instructions in SignalScan™ protocol.