<sup>ដ្</sup> ត្ត ខ្លួ រត្ត	iravir)	C T	ell Signaling
		Orders:	877-616-CELL (2355) orders@cellsignal.com
<b>寸</b> 10 mg		Support:	877-678-TECH (8324)
#5721		Web:	info@cellsignal.com cellsignal.com
#2	3 Tra	ask Lane   Danvers   Mas	sachusetts   01923   USA
	ot for Use in Diagnostic Procedures.		
Background	T-705 (Favipiravir) is a potent antiviral agent that has	s shown effectiveness aga	inst influenzas A. B. and C

	T-705 (Favipiravir) is a potent antiviral agent that has shown effectiveness against influenzas A, B, and C with IC <sub>50</sub> values as follows: A (H1N1) = 0.03-0.20 $\mu$ g/ml, A (H2N2) = 0.01-0.30 $\mu$ g/ml, A (H3N2) = 0.08-0.48 $\mu$ g/ml, A (H4N2) = 0.14-0.15 $\mu$ g/ml, A (H7N2) = 0.24-1.60 $\mu$ g/ml, B = 0.04-0.09 $\mu$ g/ml, and C = 0.03-0.06 $\mu$ g/ml (1,2). T-705 (Favipiravir) becomes activated by phosphoribosylation to form T-705 (Favipiravir) becomes activated by phosphoribosylation to form T-705 (Favipiravir)-RTP that is recognized by RNA polymerase, inhibiting RNA-dependent RNA polymerase (RdRP). This activity makes T-705 (Favipiravir) an effective inhibitor of other RNA viruses, such as arena-, phlebo-, hanta-, flavi-, entero-, and alphavirus. The antiviral capability of T-705 (Favipiravir) makes it a compound of interest when studying SARS-CoV-2 infection (3).
Molecular Formula	$C_5H_4FN_3O_2$
Molecular Weight	157.1 g/mol
Purity	>99%
CAS	259793-96-9
Solubility	Soluble in DMSO at 30 mg/ml or water at 12 mg/ml with slight warming.
Storage	Store lyophilized at -20°C, desiccated. In lyophilized form, the chemical is stable for 24 months. Once in solution, store at -20°C and use within 3 months to prevent loss of potency. <i>Aliquot to avoid multiple freeze/thaw cycles.</i>
Directions for Use	T-705 (Favipiravir) is supplied as a lyophilized powder. For a 10 mM stock, reconstitute 2 mg of powder in 1.27 ml of DMSO. Working concentrations and length of treatment can vary depending on the desired effect.
Background References	1. Furuta, Y. et al. (2013) <i>Antiviral Res</i> 100, 446-54. 2. Furuta, Y. et al. (2002) <i>Antimicrob Agents Chemother</i> 46, 977-81. 3. Dong, L. et al. (2020) <i>Drug Discov Ther</i> 14, 58-60.
Trademarks and Patents	Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.
Trademarks and Patents	Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc. All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.
Trademarks and Patents Limited Uses	All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for

respect to any third party products or services used by Customer in connection with the Products.