

PathScan® Inflammation 4-Plex Array Kit

✓ 1 Kit
(96 assays)

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

New 04/08

This product is for *in vitro* research use only and is not intended for use in humans or animals.
This product is not intended for use as a therapeutic or in diagnostic procedures.

Introduction: CST's PathScan® Inflammation 4-Plex Array Kit is based upon the sandwich ELISA principle. Rather than one immobilized antibody, these PathScan® 4-Plex Array Kits contain four different antibodies in each well of a 96-well plate. This allows the researcher to measure four different target proteins from a single well simultaneously. In addition, the incorporation of multiple phospho-specific antibody pairs allows the researcher to measure the activation state of important signaling nodes. The PathScan® Inflammation 4-Plex Array Kit allows the detection of: **Phospho-p38 α MAPK (Thr180/Tyr182)**; **Phospho-Stat3 (Tyr705)**; **Phospho-p44 MAPK (Thr202/Tyr204)**; and **Total Akt1** protein.

Briefly, 4 capture antibodies with distinct target specificity have been spotted onto the bottom of each well of a microwell plate. After incubation with cell lysates, the spotted antibodies capture the target proteins. Following extensive washing, a mixture of detection antibodies is added to detect the captured target proteins. An HRP-linked secondary solution is then used to recognize the bound detection antibodies. Chemiluminescent HRP substrate is used to produce luminescent signal.

Products Included	Quantity
Microtiter well plate with four capture antibodies per well	96 well plate
Sample diluent buffer	10 ml
20X Wash buffer	50 ml
Detection antibody mix	5.5 ml
Streptavidin-HRP solution	5.5 ml
Sealing tape	2 sheets
Chemiluminescent substrate	5.5 ml
10X Cell Lysis Buffer #9803	10 ml




Figure 1: Schematic representation of the position of targets in a single well.




Figure 2: NIH/3T3 Cells were starved and stimulated with PDGF (100 ng/ml) for 5 minutes or IL-6 (100 ng/ml) for 5 minutes at 37°C. For UV treatment, exponentially growing NIH/3T3 cells were irradiated without serum starvation. Lysates were prepared and analyzed using the PathScan® Cell-Growth 4-Plex Array Kit. Luminescent signals were acquired using a Quansys Biosciences imager. Spot intensity was quantified using Q-View™ software. The bar graph shows measurements of four different target proteins. Corresponding Western blots obtained with the same cell lysates are shown at the bottom.

Specificity/Sensitivity: PathScan® Inflammation 4-plex Array Kit detects endogenous levels of all target proteins. As shown in Figure 2, stimulation of NIH/3T3 cells with PDGF promotes phosphorylation of p44 MAPK at Thr202/Tyr204 and a modest induction of Stat3 phosphorylation on Tyr705, while p38 MAP kinase remains unaffected. Stimulation with IL-6 or UV induces phosphorylation of Stat3 or p38, respectively, while not affecting measured phosphorylation levels of the other targets. Levels of total Akt1 remain essentially unchanged with all treatments.

Species Cross-Reactivity: H, M

Required Equipment:

Quansys Imager and Q-View Plus™ software (Quansys Biosciences, UT; www.quansysbio.com) are recommended for this assay kit for acquisition of the luminescent signal and quantification of spot intensity, respectively. Alternative imagers compatible with this array product are: Alpha Innotech Fluorchem HD2 or SP; Bio-Rad Versa Doc 4000 or XRS; Fujifilm LAS-3000.

PathScan® 4-Plex Array Kit co-developed by Cell Signaling Technology® (CST) and Quansys Biosciences. Antibodies provided by CST. Platform provided by Quansys Biosciences.

Applications Key: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E—ELISA
Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebra fish B—bovine
 Dg—Dog Pg—Pig Sc—S. cerevisiae All—all species expected Species enclosed in parentheses are predicted to react based on 100% sequence homology.




Figure 3: Cell extracts from unstimulated or stimulated NIH/3T3 cells were applied to wells at increasing total protein concentrations and analyzed. Luminescent signals were acquired and spot intensity quantified. The graphs show the relationship between the total protein concentration and the resulting strength of the luminescent signal for each target.

Companion Products:

- PathScan® Total Akt1 Sandwich ELISA Kit #7170
- PathScan® Phospho-p44 MAPK (Thr202/Tyr204) Sandwich ELISA Kit #7315
- PathScan® Phospho-p38α MAPK (Thr180/Tyr182) Sandwich ELISA Kit #7140
- PathScan® Phospho-Stat3 (Tyr705) Sandwich ELISA Kit #7300
- Cell Lysis Buffer (10X) #9803

Background: Cells must respond in an appropriate fashion to many complex signaling events. Extracellular signaling cues are organized into well defined signal transduction modules that control fundamental cellular behavior. Three prominent signaling modules that are among the best characterized are the p44/42 MAP kinase (ERK MAPK), p38 Map kinase and the JAK-Stat signal transduction pathways. p44/42 MAPK is activated by a wide variety of extracellular signals including growth and neurotrophic factors, cytokines, hormones and neurotransmitters. p44/42 activation occurs through phosphorylation of threonine and tyrosine at the sequence T*EY* by a dual specificity kinase called MAP kinase kinase (MEK). p38 MAPK is dually phosphorylated in response to pro-inflammatory cytokines or cellular stress. Stat3 (Signal Transducers and Activators of Transcription) is phosphorylated at Tyr705 in response to a variety of cytokines and growth factors. Upon phosphorylation, Stat3 translocates to the nucleus where it acts as a transcription factor. Phosphorylation levels of critical molecular switches such as MAPKs therefore serve as a reliable indicator of the activation state of the entire signaling module. The profiling of phosphorylation events using phospho-specific antibodies is now widely used to investigate diagnostic pathology (1,2). While profiling of protein phosphorylation events was shown to predict the progression of a tumor to a more invasive stage (3), it has been observed that the ratio between p44/42 and p38 MAPK may predict whether tumor cells will proliferate or enter a dormant state *in vivo* (4). PathScan® Inflammation 4-Plex Array Kit provides the researcher with means to profile numerous chemical compounds and obtain in-cell relative potency (5).

Background References:

- (1) Sheehan, K.M. et al. (2008) *Oncogene* 27, 323–31.
- (2) Danna, E.A. and Nolan, G.P. (2006) *Curr Opin Chem Biol* 10, 20–7.
- (3) Irish, J.M. et al. (2006) *Nat Rev Cancer* 6, 146–55.
- (4) Aguirre-Ghiso, J.A. et al. (2003) *Cancer Res* 63, 1684–95.
- (5) Gechtman Z. (2006) *American Drug Discovery* 1:1, 44–52.

PathScan® 4-Plex Array Kit Protocol

A Reagent Preparation

1. Prior to running the assay bring the 96-well plate to room temperature.
2. Dilute the 20X Wash Buffer to a 1X solution with deionized water. Mix thoroughly. Store at 4°C until ready to use.
3. Thaw the 1X Cell Lysis Buffer (#9803) and mix gently but thoroughly. Supplement the cell lysis buffer with PMSF to 1 mM final concentration prior to performing the assay (other protease inhibitors may be added).

B Preparing Cell Lysates

1. Place cells on ice and wash once with ice-cold PBS.
2. Scrape cells in 0.5 ml of ice-cold cell lysis buffer (per 10–15 cm plate), gently mix and incubate for 5 minutes on ice.
3. Sonicate cells on ice using medium strength output by creating 4 sonic pulses, 10 seconds each, with 30–60 sec. intervals.
4. Microcentrifuge for 10 minutes at 4°C to pellet the insoluble material and transfer the supernatant to a new tube. The resultant supernatant is the cell lysate that is ready to be analyzed. The cell lysate can also be stored at -80°C in single-use aliquots for future use.
5. Determine protein concentration of the cell lysates and make appropriate dilutions using the sample diluent buffer. Recommended range of final total protein concentration to be applied onto 96-well plates is: 0.1–0.5 mg/ml. (Use the dose response curves that appear in the datasheet as a reference). As a guideline, the total protein concentration of a typical cell lysate generated from adherent cells grown on a 15 cm plate to 80–90 % confluency will be in the range of 0.5–2.5 mg/ml.

C Test Procedure

1. Add 50 µl of diluted cell extract into each well.
2. Cover the 96-well plate with plate seal and incubate on a shaker (~120 rpm) for 1 hour at room temperature.
3. Wash the plate 3 times with Wash Buffer, 200–400 µl per well.
4. Add 50 µl of the Detection mix into each well.
5. Reseal the plate with the provided seal and incubate on a shaker (~120 rpm) for 1 hour at room temperature.
6. Wash the plate 3 times with Wash Buffer.
7. Add 50 µl of the Streptavidin-HRP solution into each well.
8. Reseal the plate with a plate seal and incubate on a shaker (~120 rpm) for 15 min. at room temperature.
9. Prepare the chemiluminescent substrate by mixing the Substrate bottle A with Substrate bottle B and set aside at room temperature.
10. Wash the plate 6 times with Wash Buffer.
11. Add 50 µl of the mixed chemiluminescent substrate into each well.
12. Acquire Luminescent signal by imaging the plate using a Quansys imager (Quansys Biosciences, UT; www.quansysbio.com; Cat # 104450GR; www.quansysbio.com/support/imagers/quansys.html). Spot intensity can be quantified using Q-View™ or Q-View Plus™ softwares from Quansys Biosciences. Alternative imagers compatible with this array product are: Alpha Innotech Fluorchem HD2 or SP; Bio-Rad Versa Doc 4000 or XRS; or Fujifilm LAS-3000.