

# GSK-3 $\beta$ Kinase

✓ 5  $\mu$ g



Cell Signaling  
TECHNOLOGY®

**Orders** ■ 877-616-CELL (2355)  
orders@cellsignal.com  
**Support** ■ 877-678-TECH (8324)  
info@cellsignal.com  
**Web** ■ www.cellsignal.com

new 07/06

This product is for *in vitro* research use only and is not intended for use in humans or animals.

**Description:** Purified recombinant full length human GSK-3 $\beta$  (Met1-Thr433) kinase, supplied as a GST fusion protein.

**Background:** Glycogen synthase kinase-3 (GSK-3) was initially identified as an enzyme that regulates glycogen synthesis in response to insulin (1). GSK-3 is a ubiquitously expressed serine/threonine protein kinase that phosphorylates and inactivates glycogen synthase. GSK-3 is a critical downstream element of the PI3 kinase/Akt cell survival pathway, and its activity can be inhibited by Akt-mediated phosphorylation at Ser21 of GSK-3 $\alpha$  and Ser9 of GSK-3 $\beta$  (2,3). GSK-3 has been implicated in the regulation of cell fate in Dictyostelium, and is a component of the Wnt signaling pathway required for *Drosophila*, *Xenopus* and mammalian development (4). GSK-3 has been shown to regulate cyclin D1 proteolysis and subcellular localization (5).

**Source/Purification:** The GST-Kinase fusion protein was produced using a baculovirus expression system with a construct expressing human GSK-3 $\beta$  (Met1-Thr433) (GenBank Accession No. NM\_002093) with an amino-terminal GST tag. The protein was purified by one-step affinity chromatography using glutathione-agarose.

**Quality Control:** The theoretical molecular weight of the GST-GSK-3 $\beta$  fusion protein is 73 kDa. The purified kinase was quality controlled for purity using SDS-PAGE followed by Coomassie stain [Fig.1]. GSK-3 $\beta$  kinase activity was determined using a radiometric assay [Fig.2].




Figure 1. The purity of the GST-GSK-3 $\beta$  fusion protein was analyzed using SDS/PAGE followed by Coomassie stain.

## Background References:

- (1) Welsh, G.I. et al. (1996) *Trends Cell Biol.* 6, 274–279.
- (2) Srivastava, A.K. and Pandey, S.K. (1998) *Mol. Cell. Biochem.* 182, 135–141.
- (3) Cross, D.A. et al. (1995) *Nature* 378, 785–789.
- (4) Nusse, R. (1997) *Cell* 89, 321–323.
- (5) Diehl, J.A. et al. (1998) *Genes Dev.* 12, 3499–3511.




Figure 2. GSK-3 $\beta$  kinase activity was measured in a radiometric assay using the following reaction conditions: 4 mM MOPS, pH 7.2, 2.5 mM  $\beta$ -glycerophosphate, 1 mM EGTA, 0.4 mM EDTA, 4 mM MgCl<sub>2</sub>, 0.05 mM DTT, 40 ng/ $\mu$ l BSA, 50  $\mu$ M ATP, Substrate: GSK3 substrate peptide 400 ng/ $\mu$ l and recombinant GSK-3 $\beta$ : variable.

**Storage:** Enzyme is supplied in 50 mM Tris-HCl, pH 7.5; 150 mM NaCl, 0.25 mM DTT, 0.1 mM EGTA, 0.1 mM EDTA, 0.1 mM PMSF, 25% glycerol, 7 mM glutathione. Store at -80°C.

Keep on ice during use.

Avoid repeated freeze-thaw cycles.

## Protocol for GSK-3 $\beta$ Kinase Assay

**Note:** Lot-specific information for this kinase is provided on the enzyme vial. Optimal assay incubation times and enzyme concentrations must be determined empirically for each lot of kinase under specified conditions.

### A Additional Solutions and Reagents (Not included)

#### 1. Kinase Buffer (10X)

50 mM MOPS, pH 7.2  
25 mM  $\beta$ -glycerophosphate  
10 mM EGTA  
4 mM EDTA  
40 mM  $\text{MgCl}_2$   
0.5 mM DTT  
400 ng/ $\mu$ l BSA

#### 2. ATP (10 mM) #9804

#### 3. $^{32}\text{P}$ - $\gamma$ ATP

#### 4. GSK substrate (YRRAAVPPSPSLSRHSSPHQS\*EDEEE) (1 mg/ml)

### B Suggested Protocol

1. Dilute 10 mM ATP with 3X assay buffer 1:40 to make 250  $\mu$ M ATP.
2. Dilute [ $^{32}\text{P}$ ] ATP to 0.16  $\mu\text{Ci}/\mu\text{l}$  [ $^{32}\text{P}$ ] ATP with 250  $\mu$ M ATP solution.
3. Transfer enzyme from  $-80^\circ\text{C}$  to ice. Allow enzyme to thaw on ice.
4. Dilute GSK-3 $\beta$  protein (100 ng/ $\mu$ l concentration) to 20 ng/ $\mu$ l with 1X assay buffer followed by 2-fold serial dilutions.
5. To start the reaction combine 10  $\mu$ l diluted GSK-3 $\beta$  kinase solution, 10  $\mu$ l GSK3 substrate peptide (1  $\mu\text{g}/\mu\text{l}$ ) and 5  $\mu$ l 0.16  $\mu\text{Ci}/\mu\text{l}$  [ $^{32}\text{P}$ ] ATP solution.

### Final Assay Conditions

4 mM MOPS, pH 7.2  
2.5 mM  $\beta$ -glycerophosphate  
1 mM EGTA  
0.4 mM EDTA  
4 mM  $\text{MgCl}_2$   
0.05 mM DTT  
40  $\mu$ M BSA  
400 ng/ $\mu$ l GSK substrate

6. After 15 minutes terminate reaction by spotting 20  $\mu$ l of the reaction mixture onto phosphocellulose P81 paper.
7. Air dry the P81 paper then wash with 1% phosphoric acid 3 times.
8. Transfer P81 paper to 4 ml scintillation tube then add 3 ml scintillation cocktail.
9. Count samples in a scintillation counter.

Cell Signaling Technology offers a full line of protein kinases, substrates, and antibody detection reagents for high throughput screening. Please direct all inquiries to: [drugdiscovery@cellsignal.com](mailto:drugdiscovery@cellsignal.com).