# HTScan® Pim-1 Kinase Assay Kit

100 assays (96 Well Format)



OrdersImage: 877-616-CELL (2355)<br/>orders@cellsignal.comSupportImage: 877-678-TECH (8324)<br/>info@cellsignal.comWebImage: www.cellsignal.com

rev. 08/13/07

This product is for in vitro research use only and is not intended for use in humans or animals.

| Products Included                 | Products # | Kit Quantity |
|-----------------------------------|------------|--------------|
| Phospho-Bad (Ser112) Antibody     | 9291       | 30 µl        |
| Kinase Buffer (10X)               | 9802       | 15 ml        |
| ATP (10 mM)                       | 9804       | 1 ml         |
| Bad (Ser112) Biotinylated Peptide | 1342       | 1.25 ml      |
| Pim-1 Kinase (recombinant, human) | 7572       | 5 µg         |

**Description:** The kit provides a means of performing kinase activity assays with recombinant human Pim-1 kinase. It includes active Pim-1 kinase (supplied as a GST fusion protein), a biotinylated peptide substrate and a phospho-serine/threonine antibody for detection of the phosphorylated form of the substrate peptide.

#### Peptide Core Sequence: RS\*RHS\*S\*Y

**Molecular Weights:** Peptide substrate, Biotin-Bad (Ser112): 1,700 Daltons. GST-Pim-1 Kinase: 65 kDa.

**Source/Purification:** The GST-Kinase fusion protein was produced using a baculovirus expression system with a construct expressing full length human Pim-1 (Met1-Lys313) (GenBank Accession No. NM\_002648) with an amino-terminal GST tag. The protein was purified by onestep affinity chromatography using glutathione-agarose. Quality Control: The substrate peptide was selected using our Serine/Threonine Kinase Substrate Screening Kit #7400. Phospho-Bad (Ser112) Antibody #9291 was used for detection. The quality of the biotinylated peptide was evaluated by reverse-phase HPLC and by mass spectrometry.

Purified Pim-1 kinase was quality controlled for purity by SDS-PAGE followed by Coomassie stain and Western blot. The specific activity of the Pim-1 was determined using a radiometric assay [Fig.1]. Time course [Fig.2], kinase dose dependency [Fig.3] and substrate dose-dependency [Fig.4] assays were performed to verify Pim-1 activity using the Pim-1 substrate peptide provided in this kit. Pim-1 sensitivity to the inhibitor staurosporine was measured using the Pim-1 substrate peptide provided in this kit [Fig.5].

**Storage:** Antibodies are supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100  $\mu$ g/ml BSA and 50% glycerol. Do not aliquot the antibodies. Peptides are supplied at 6  $\mu$ M in 0.001% DMSO. Enzymes are supplied in 50 mM Tris-HCL (pH 8.0), 100 mM NaCl, 5 mM DTT, 15 mM reduced glutathione and 20% glycerol. Store at  $-80^{\circ}$ C.

Keep enzymes on ice during use.

Avoid repeated freeze-thaw cycles.

#### **Companion Products:**

Serine/Threonine Kinase Substrate Screening Kit #7400

Pim-1 Kinase #7572

Phospho-Bad (Ser112) Antibody #9291

Bad (Ser112) Biotinylated Peptide #1342

Staurosporine #9953

Kinase Buffer (10X) #9802

ATP (10 mM) #9804




Figure 1. Pim-1 kinase activity was measured in a radiometric assay using the following reaction conditions: 60 mM HEPES-NaOH, pH 7.5, 3 mM MgCl<sub>2</sub>, 3 mM MnCl<sub>2</sub>, 3 µM Na-orthovanadate, 1.2 mM DTT, ATP (variable), 2.5 µg/50 µl PEG20,000, Substrate: Histone H2B 5 µg/50 µl, recombinant Pim-1: 200 ng/50 µl.






Figure 2. Time course of Pim-1 kinase activity: DELFIA® data generated using Phospho-Bad (Ser112) Antibody #9291 to detect phosphorylation of Pim-1 substrate peptide (#1342) by Pim-1 kinase. In a 50 µl reaction, 50 ng Pim-1 and 1.5 µM substrate peptide were used per reaction. (DELFIA® is a registered trademark of PerkinElmer, Inc.)




Figure 3. Dose dependence curve of Pim-1 kinase activity: DELFIA® data generated using Phospho-Bad (Ser112) Antibody #9291 to detect phosphorylation of substrate peptide (#1342) by Pim-1 kinase. In a 50 µl reaction, increasing amounts of Pim-1 and 1.5 µM substrate peptide were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)




Figure 4. Peptide concentration dependence of Pim-1 kinase activity: DELFIA® data generated using Phospho-Bad (Ser112) Antibody #9291 to detect phosphorylation of substrate peptide (#1342) by Pim-1 kinase. In a 50 µl reaction, 50 ng of Pim-1 and increasing concentrations of substrate peptide were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)




Figure 5. Staurosporine inhibition of Pim-1 kinase activity: DELFIA® data generated using Phospho-Bad (Ser112) Antibody #9291 to detect phosphorylation of Pim-1 substrate peptide (#1342) by Pim-1 kinase. In a 50 µl reaction, 50 ng Pim-1, 1.5 µM substrate peptide, 20 µM ATP and increasing amounts of staurosporine were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)

Background: Pim proteins (Pim-1, Pim-2 and Pim-3) are oncogene-encoded serine/threonine kinases (1). Pim-1, a serine/threonine kinase highly expressed in hematopoietic cells, plays a critical role in the transduction of mitogenic signals and is rapidly induced by a variety of growth factors and cytokines (1-4). Pim-1 cooperates with c-Myc in lymphoid cell transformation and protects from growth factor withdrawal and genotoxic stress-induced apoptosis (5,6). Pim-1 also enhances the transcriptional activity of c-Myb through direct phosphorylation within the c-Myb DNA binding domain as well as phosphorylation of the transcriptional coactivator p100 (7,8). Hypermutations of the Pim-1 gene are found in B-cell diffuse large cell lymphomas (9). Phosphorylation of Pim-1 at Tyr 218 by Etk occurs following IL-6 stimulation and is correlated with an increase in Pim-1 activity (10). Various substrates have been identified for Pims. Among them, BAD has been shown to be phosphorylated by both Pim-1 and Pim-2 at Ser112 and this phosphorylation reverses BAD-induced cell apoptosis (11,12).

#### **Background References:**

(1) Mikkers, H. et al. (2004) *Mol. Cell.Biol.* 24, 6104–6115.

Cell Signaling

TECHNOLOGY®

- (2) Selten, G. et al. (1986) Cell 46, 603-611.
- (3) Meeker, T.C. et al. (1987) J. Cell. Biochem. 35, 105–112.
- (4) Dautry, F. et al. (1988) *J. Biol. Chem.* 263, 17615–17620.
- (5) Moroy, T. et al. (1993) Proc. Natl. Acad. Sci. USA 90, 10734–10738.
- (6) Lilly, M. and Kraft, A. (1997) *Cancer Res.* 57, 5348–5355.
- (7) Leverson, J.D. et al. (1998) Mol. Cell. 2, 417-425.
- (8) Winn, L.M. et al. (2003) Cell Cycle 2, 258-262.
- (9) Pasqualucci, L. et al. (2001) Nature 412, 341-346.
- (10) Kim, O. et al. (2004) Oncogene 23, 1838-1844.
- (11) Aho, T. L. et al. (2004) FEBS Lett. 571, 43-49.
- (12) Yan, B. et al. (2003) *J. Biol. Chem.* 278, 45358–45367.

# Protocol for HTScan® Pim-1 Kinase Assay Kit

#### Kinase

**Note:** Lot-specific information for this kinase is provided on the enzyme vial. Optimal assay incubation times and enzyme concentrations must be determined empirically for each lot of kinase under specified conditions.

## A Additional Solutions and Reagents (Not included)

- 1. Wash Buffer: 1X PBS, 0.05% Tween-20 (PBS/T)
- **2.** Bovine Serum Albumin (BSA)
- 3. Stop Buffer: 50 mM EDTA pH 8

DELFIA® is a registered trademark of PerkinElmer Life Sciences

#### B Suggested Protocol for 100 Assays

- Add 100 µl 10 mM ATP to 1.25 ml 6 µM substrate peptide. Dilute the mixture with dH<sub>2</sub>0 to 2.5 ml to make 2X ATP/substrate cocktail ([ATP]=400 µM, [substrate] = 3 µm).
- 2. Transfer enzyme from -80°C to ice. Allow enzyme to thaw on ice.
- Microcentrifuge briefly at 4°C to bring liquid to the bottom of the vial. Return immediately to ice.
- Add 1 ml 10X kinase buffer [1 ml 10X Kinase Buffer 250 mM Tris-HCl pH 7.5, 100 mM MgCl<sub>2</sub>, 1 mM Na<sub>3</sub>VO<sub>4</sub>, 50 mM β-glycerophosphate, 20 mM dithiothreitol (DTT)] to 1.5 ml dH<sub>2</sub>O to make 2.5 ml 4X reaction buffer.
- Transfer 1.2 ml of 4X Reaction but ffer to each enzyme tube to make 4X reaction cocktail ([enzyme]) = 4 ng/µl in 4X reaction cocktail).
- Add 12.5 µl of the 4X reaction cocktail to 12.5 µl/well of prediluted compound of interest (usually around 10 µM) and incubate for 5 minutes at room temperature.
- **7.** Add 25 μl of 2X ATP/substrate cocktail to 25 μl/well preincubated reaction cocktail/compound.

#### Final Assay Conditions for a 50 $\mu I$ Reaction

 $\begin{array}{l} \text{25 mM Tris-HCI (pH 7.5)} \\ \text{10 mM MgCl}_2 \\ \text{5 mM } \beta \text{-glycerophosphate} \\ \text{0.1 mM Na}_3 \text{VO}_4 \\ \text{2 mM DTT} \\ \text{200 } \mu\text{M ATP} \\ \text{1.5 } \mu\text{M peptide} \\ \text{50 ng Pim-1 Kinase} \end{array}$ 

- 8. Incubate reaction plate at room temperature for 30 minutes.
- **9.** Add 50  $\mu$ /well Stop Buffer (50 mM EDTA, pH 8) to stop the reaction.
- **10.** Transfer 25  $\mu$ l of each reaction to a 96-well streptavidin-coated plate containing 75  $\mu$ l dH<sub>2</sub>O/well and incubate at room temperature for 60 minutes.
- **11.** \*Wash three times with 200  $\mu$ I/well PBS/T.
- Dilute primary antibody, Phospho-Bad (Ser112) Antibody #9291, 1:1000 in PBS/T with 1% BSA. Add 100 μl/well primary antibody.
- **13.** Incubate at room temperature for 120 minutes.
- **14.** \*Wash three times with 200 µl/well PBS/T.
- 15. For  $\mathsf{DELFIA}^{\otimes}$  or Colorimetric ELISA detection methods please use the following protocols.

## **DELFIA®** Assay

- Prepare appropriate dilution of Europium labeled secondary antibody in PBS/T with 1% BSA (1:500 dilution for anti-mouse IgG or 1:1000 for anti-rabbit IgG).
- 2. Add 100  $\mu\text{I/well}$  secondary antibody solution.
- **3.** Incubate at room temperature for 30 minutes.
- 4. \*Wash five times with 200  $\mu$ l/well PBS/T.
- 5. Add 100  $\mu\text{I/well}$  DELFIA® Enhancement Solution.
- 6. Incubate at room temperature for 5 minutes.
- Read plate using a Time Resolved Fluorescent plate reader using the following settings;
  - a. Excitation Filter: 340 nm
  - b. Emission Filter: 615 nm
  - c. Delay\*\*: 400 µs
- \*\* Delay time is the delay from the excitation pulse to the beginning of the measurement.

#### **Companion Products for DELFIA®**

DELFIA® Europium-labeled Anti-mouse IgG (PerkinElmer Life Sciences #AD0124) DELFIA® Europium-labeled Anti-rabbit IgG (PerkinElmer Life Sciences #AD0105) DELFIA® Enhancement Solution (PerkinElmer Life Sciences #1244-105) DELFIA® Streptavidin coated, 96-well, yellow plate (PerkinElmer Life Sciences AAAND-0005)

# **Colorimetric ELISA Assay**

- Prepare appropriate dilution of HRP labeled secondary antibody in PBS/T with 1% BSA (1:500 dilution for anti-mouse IgG or 1:1000 for anti-rabbit IgG).
- 2. Add 100 µl/well secondary antibody solution.
- **3.** Incubate at room temperature for 30 minutes.
- 4. \*Wash five times with 200 µl/well PBS/T.
- 5. Add 100 µl/well TMB substrate.
- **6.** Incubate at room temperature for 15 minutes.
- 7. Add 100 µl/well of stop solution.
- 8. Mix well.
- 9. Read the absorbance at 405 nm with a microtiter plate reader.

#### **Companion Products For Colorimetric ELISA Assay**

Anti-mouse IgG, HRP Linked Antibody #7076 Anti-rabbit IgG, HRP Linked Antibody #7074 TMB Solution #7004 Stop Solution #7002

\*NOTE: Use of an automated microplate washer as well as centrifugation of plates when appropriate, greatly improves reproducibility.

Please contact Cell Signaling Technology for HTS-ready antibodies (PBS formulated and carrier-free), and detailed peptide substrate sequence information. Email: drugdiscovery@cellsignal.com