Jak2 **Kinase**

☑ 5 µg



Orders 877-616-CELL (2355)

orders@cellsignal.com

Support 877-678-TECH (8324)

info@cellsignal.com

Web www.cellsignal.com

rev. 12/06/05

This product is for in vitro research use only and is not intended for use in humans or animals.

Description: Purified recombinant murine Jak2 kinase (Lys752-Val1129), supplied as a GST fusion protein.

Background: The Janus family of tyrosine kinases, consisting of Jak1, Jak2, Jak3 and Tyk2, associates with a variety of cytokine receptors and are activated by ligand binding to these receptors (1). Upon ligand binding to cytokine receptors, Jaks autophosphorylate and phosphorylate their associated receptors, providing multiple binding sites for signaling proteins containing SH2 or other phosphotyrosine-binding domains, including Stats (2), Shc (3), insulin receptor substrates (4) and focal adhesion kinase (FAK) (5).

Source/Purification: The GST-Jak2 fusion protein was produced using a baculovirus expression system with a construct expressing murine Jak2 (Lys752-Val1129) (GenBank accession No. L16956) with an amino-terminal GST tag. The protein was purified by one-step affinity chromatography using glutathione-agarose.

Quality Control: The theoretical molecular weight of the GST-Jak2 fusion protein is 71 kDa. The purified kinase was quality controlled for purity using SDS-PAGE followed by Silver stain and Western blot [Fig.1]. Jak2 kinase activity was determined using a radiometric assay [Fig.2]. A kinase dose dependency assay was performed to measure Jak2 activity using HTScan® Jak2 Kinase Assay Kit #7752 [Fig.3].

Background References:

- (1) Leonard, W.J. and O'Shea, J.J. (1998) Annu. Rev. Immunol. 16, 293-322.
- (2) Darnell, J.E. (1997) Science 277, 1630-5.
- (3) VanderKuur, J. et al. (1995) J. Biol. Chem. 270, 7587-93.
- (4) Argetsinger, L.S. et al. (1995) J. Biol. Chem. 270, 14685-92.
- (5) Zhu, T. et al. (1998) J. Biol. Chem. 273, 10682-9.

Storage: Enzyme is supplied in 50 mM Tris-HCl, pH 8.0; 100 mM NaCl, 5 mM DTT, 15 mM reduced glutathione, 20% glycerol. Store at -80°C.

Keep on ice during use.

Avoid repeated freeze-thaw cycles.

Companion Products:

HTScan® Jak2 Kinase Assay Kit #7752

FLT3 (Tyr589) Biotinylated Peptide #1305

Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411

Kinase Buffer (4X) #9805

ATP (10 mM) #9804

Staurosporine #9953

Tyrosine Kinase Substrate Screening Kit #7450




Figure 1. The purity of the GST-Jak2 fusion protein was analyzed using SDS/PAGE followed by anti-Jak2 Western blot (A) or Silver stain (B).




Figure 2. Jak2 kinase activity was measured in a radiometric assay using the following reaction conditions: 60 mM HEPES-NaOH, pH 7.5, 3 mM MgCl₂, 3 mM MnCl₂, 3 µM Na-orthovanadate, 1.2 mM DTT, ATP (variable), 2.5 µg/50 µl PEG20,000, Substrate: PolyEY, 1.5 µg/50 µl and 50 ng/50 µl Recombinant Jak2.




Figure 3. Dose dependence curve of Jak2 kinase activity: DELFIA® data generated using Phospho-Tyrosine mAb #9411 to detect phosphorylation of substrate peptide (#1305) by Jak2 kinase. In a 50 µI reaction, increasing amounts of Jak2 and 1.5 µM substrate peptide were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)



Protocol for Jak2 Kinase Assay

*IMPORTANT: Use of an automated microplate washer as well as centrifugation of plates when appropriate, greatly improves reproducibility.

Kinase

Note: Lot-specific information for this kinase is provided on the enzyme vial. Optimal assay incubation times and enzyme concentrations must be determined empirically for each lot of kinase under specified conditions.

Additional Solutions and Reagents (Not included)

- 1. Wash Buffer: 1X PBS, 0.05% Tween-20 (PBS/T)
- 2. Bovine Serum Albumin (BSA)
- Stop Buffer: 50 mM EDTA pH 8
- 4. Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411
- 5. HTScan® Tyrosine Kinase Buffer (4X) #9805
- 6. ATP (10 mM) #9804
- 7. FLT3 (Tyr589) Biotinylated Peptide #1305
- 8. DTT (1000X, 1.25 M)
- 9. DELFIA® Europium-labeled Anti-mouse IgG (PerkinElmer Life Sciences
- 10. DELFIA® Enhancement Solution (PerkinElmer Life Sciences #1244-105)
- 11. DELFIA® Streptavidin coated, 96-well, yellow plate (PerkinElmer Life Sciences AAAND-0005)

DELFIA® is a registered trademark of PerkinElmer Life Sciences

Suggested Protocol for 100 Assays

- 1. Add 100 μ l 10 mM ATP to 1.25 ml 6 μ M substrate peptide. Dilute the mixture with dH₂0 to 2.5 ml to make 2X ATP/substrate cocktail ([ATP]=400 μM, [substrate] = $3 \mu m$).
- 2. Immediately transfer enzyme from -80°C to ice. Allow enzyme to thaw on ice.
- 3. Microcentrifuge briefly at 4°C to bring liquid to the bottom of the vial. Return immediately to ice.
- 4. Add 10 µl of DTT (1.25 M) to 2.5 ml of 4X HTScan™ Tyrosine Kinase Buffer (240 mM HEPES pH 7.5, 20 mM MgCl₂, 20 mM MnCl₂, 12 μM Na₂VO₄) to make DTT/Kinase buffer.
- **5.** Transfer 1.25 ml of DTT/Kinase buffer to enzyme tube to make 4X reaction cocktail ([enzyme]=8.0 ng/µl in 4X reaction cocktail).
- **6.** Incubate 12.5 μl of the 4X reaction cocktail with 12.5 μl/well of prediluted compound of interest (usually around 10 µM) for 5 minutes at room temperature.
- 7. Add 25 µl of 2X ATP/substrate cocktail to 25 µl/well preincubated reaction cocktail/compound.

Final Assay Conditions for a 50 µl Reaction

60 mM HEPES pH 7.5

5 mM MgCl₂

5 mM MnCl

3 µM Na₂VO₄

1.25 mM DTT

200 uM ATP 1.5 µM peptide

100 ng Jak2 Kinase

- 8. Incubate reaction plate at room temperature for 30 minutes.
- 9. Add 50 µl/well Stop Buffer (50 mM EDTA, pH 8) to stop the reaction.
- 10. Transfer 25 µl of each reaction and 75 µl of dH₂0/well to a 96-well streptavidincoated plate and incubate at room temperature for 60 minutes.
- 11. *Wash three times with 200 µl/well PBS/T.
- 12. Dilute primary antibody in PBS/T with 1% BSA. Add 100 µl/well primary antibody.

Please note: This protocol was validated using a FLT3 (Tyr589) Biotinylated Peptide and Phospho-Tyrosine Mouse mAb (P-Tyr-100) diluted 1:1000 (see additional reagents). Primary antibody chosen should be specific to the

- **13.** Incubate at room temperature for 120 minutes.
- 14. *Wash three times with 200 µl/well PBS/T.
- 15. Dilute Europium labeled anti-mouse IgG 1:500 in PBS/T with 1% BSA. Add 100 µl/well diluted antibody.
- **16.** Incubate at room temperature for 30 minutes.
- 17. *Wash five times with 200 µl/well PBS/T.
- 18. Add 100 µl/well DELFIA® Enhancement Solution.
- 19. Incubate at room temperature for 5 minutes.
- 20. Detect 615 nm fluorescence emission with appropriate Time-Resolved Plate Reader.

Please contact Cell Signaling Technology for HTS-ready antibodies (PBS formulated and carrier-free), and detailed peptide substrate sequence information.

Email: drugdiscovery@cellsignal.com