# HTScan® VEGF Receptor 3 Kinase Assay Kit

100 assays (96 Well Format)



Orders=877-616-CELL (2355)<br/>orders@cellsignal.comSupport=877-678-TECH (8324)<br/>info@cellsignal.comWeb=www.cellsignal.com

rev. 09/05/07

This product is for in vitro research use only and is not intended for use in humans or animals.

Products Included	Products #	Kit Quantity
Phospho-Tyrosine Mouse mAb (P-Tyr-100)	9411	30 µl
HTScan® Tyrosine Kinase Buffer (4X)	9805	15 ml
DTT (1000X)		80 µl
ATP (10 mM)	9804	1 ml
MET (Tyr1253) Biotinylated Peptide	1367	1.25 ml
VEGF Receptor 3 Kinase (recombinant, human)	7790	2 X 5 µg

**Description:** The kit provides a means of performing kinase activity assays with recombinant human VEGFR-3 kinase. It includes active VEGFR-3 kinase (supplied as a GST fusion protein), a biotinylated peptide substrate and a phospho-tyrosine antibody for detection of the phosphorylated form of the substrate peptide.

## Peptide Core Sequence: MY\*DKEY\*Y\*S

**Molecular Weights:** Peptide substrate, Biotin-MET (Tyr1253) peptide: 2,200Daltons. GST-VEGFR-3 Kinase: 86 kDa.

**Background:** Vascular endothelial growth factor receptor 3 (VEGFR-3) is a 195 kDa membrane receptor tyrosine kinase. All of the VEGF receptors are characterized by seven extracellular immunoglobulin (Ig)-like domains followed by a membrane-spanning domain and a conserved intracellular tyrosine kinase domain (1). VEGFR-3 is largely restricted to the lymphatic endothelium in adult tissue and is thought to control lymphangiogenesis (1,2). Binding of VEGF-C/VEGF-D to VEGFR-3 results in transphosphorylation of tyrosine residues in its intracellular domain, recruitment of signaling molecules and activation of ERK1/2 and Akt signaling cascades (1,3).



Figure 1. VEGFR-3 kinase activity was measured in a radiometric assay using the following reaction conditions: 60 mM HEPES-NaOH, pH 7.5, 3 mM MgCl<sub>2</sub>, 3 mM MnCl<sub>2</sub>, 3 µM Na-orthovanadate, 1.2 mM DTT, ATP (variable), 2.5 µg/50 µl PEG20,000, Substrate: PolyEY5 µg/50 µl, recombinant VEGFR-3: 100 ng/50 µl.

**Source/Purification:** The GST-Kinase fusion protein was produced using a baculovirus expression system with a construct expressing a fragment of human VEGFR-3 (Asn799-Arg1298) (GenBank accession No. NM\_002020) with an amino-terminal GST tag. The protein was purified by one-step affinity chromatography using glutathioneagarose.

Quality Control: The substrate peptide was selected using our Tyrosine Kinase Substrate Screening Kit #7450. Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411 was used for detection. The quality of the biotinylated peptide was evaluated by reverse-phase HPLC and by mass spectrometry.

Purified VEGFR-3 kinase was quality controlled for purity by SDS-PAGE followed by Coomassie stain and Western blot. The specific activity of the VEGFR-3 kinase was determined using a radiometric assay [Fig.1]. Time course [Fig.2], kinase dose dependency [Fig.3] and substrate dose-dependency [Fig.4] assays were performed to verify VEGFR-3 activity using the VEGFR-3 substrate peptide provided in this kit. VEGFR-3 sensitivity to the inhibitor staurosporine was measured using the VEGFR-3 substrate peptide provided in this kit [Fig.5].

# Background References:

(1) Robinson, C.J. and Stringer, S.E. (2001) *J. Cell Sci.*114, 853–865.

- (2) Valtola, R. et al. (1999) Am. J. Pathol. 154, 1381–1390.
- (3) Saharinen, P. and Petrova, T.V. (2004) Ann. N.Y. Acad. Sci. 1014, 76–87.

**Storage:** Antibodies are supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100  $\mu$ g/ml BSA and 50% glycerol. Do not aliquot the antibodies. Peptides are supplied at 6  $\mu$ M in 0.001% DMSO. Enzymes are supplied in 50 mM Tris-HCL (pH 8.0), 100 mM NaCl, 5 mM DTT, 15 mM reduced glutathione and 20% glycerol. Store at  $-80^{\circ}$ C.

Keep enzymes on ice during use.

Avoid repeated freeze-thaw cycles.

### **Companion Products:**

Tyrosine Kinase Substrate Screening Kit #7450

VEGF Receptor 3 Kinase #7790

Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411

MET (Tyr1253) Biotinylated Peptide #1367

Staurosporine #9953







Figure 2. Time course of VEGFR-3 kinase activity: DELFIA® data generated using Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411 to detect phosphorylation of VEGFR-3 substrate peptide (#1367) by VEGFR-3 kinase. In a 50 µl reaction, 100 ng VEGFR-3 and 1.5 µM substrate peptide were used per reaction. (DELFIA® is a registered trademark of PerkinElmer, Inc.)



Figure 3. Dose dependence curve of VEGFR-3 kinase activity: DELFIA® data generated using Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411 to detect phosphorylation of substrate peptide (#1367) by VEGFR-3 kinase. In a 50 µl reaction, increasing amounts of VEGFR-3 and 1.5 µM substrate peptide were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)



Figure 4. Peptide concentration dependence of VEGFR-3 kinase activity: DELFIA® data generated using Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411 to detect phosphorylation of substrate peptide (#1367) by VEGFR-3 kinase. In a 50 µl reaction, 100 ng of VEGFR-3 and increasing concentrations of substrate peptide were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)



Figure 5. Staurosporine inhibition of VEGFR-3 kinase activity: DELFIA® data generated using Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411 to detect phosphorylation of VEGFR-3 substrate peptide (#1367) by VEGFR-3 kinase. In a 50 µl reaction, 100 ng VEGFR-3, 1.5 µM substrate peptide, 20 µM ATP and increasing amounts of staurosporine were used per reaction at room temperature for 30 minutes. (DELFIA® is a registered trademark of PerkinElmer, Inc.)

# Protocol for HTScan® VEGF Receptor 3 Kinase Assay Kit

# Kinase

**Note:** Lot-specific information for this kinase is provided on the enzyme vial. Optimal assay incubation times and enzyme concentrations must be determined empirically for each lot of kinase under specified conditions.

# A Additional Solutions and Reagents (Not included)

- 1. Wash Buffer: 1X PBS, 0.05% Tween-20 (PBS/T)
- **2.** Bovine Serum Albumin (BSA)
- 3. Stop Buffer: 50 mM EDTA pH 8

DELFIA® is a registered trademark of PerkinElmer Life Sciences

# B Suggested Protocol for 100 Assays

- Add 10 µl 10 mM ATP to 1.25 ml 6 µM substrate peptide. Dilute the mixture with dH<sub>2</sub>0 to 2.5 ml to make 2X ATP/substrate cocktail ([ATP]=40 µM, [substrate]=3 µm).
- 2. Immediately transfer enzyme from -80°C to ice. Allow enzyme to thaw on ice.
- Microcentrifuge briefly at 4°C to bring liquid to the bottom of the vial. Return immediately to ice.
- Add 10 μl of DTT (1.25 M) to 2.5 ml of 4X HTScan<sup>®</sup> Tyrosine Kinase Buffer (240 mM HEPES pH 7.5, 20 mM MgCl<sub>2</sub>, 20 mM MnCl<sub>2</sub>, 12 μM Na<sub>3</sub>VO<sub>4</sub>) to make DTT/Kinase buffer.
- Transfer 0.6 ml of DTT/Kinase buffer to each enzyme tube to make 4X reaction cocktail ([enzyme] = 8 ng/µL in 4X reaction cocktail).
- Incubate 12.5 μl of the 4X reaction cocktail with 12.5 μl/well of prediluted compound of interest (usually around 10 μM) for 5 minutes at room temperature.
- Add 25 μl of 2X ATP/substrate cocktail to 25 μl/well preincubated reaction cocktail/compound.

#### Final Assay Conditions for a 50 $\mu I$ Reaction

 $\begin{array}{l} \text{60 mM HEPES pH 7.5} \\ \text{5 mM MgCl}_2 \\ \text{5 mM MnCl}_2 \\ \text{3 } \mu\text{M Na}_3\text{VO}_4 \\ \text{1.25 mM DTT} \\ \text{20 } \mu\text{M ATP} \\ \text{1.5 } \mu\text{M peptide} \\ \text{100 ng VEGF Receptor 3 Kinase} \end{array}$ 

- 8. Incubate reaction plate at room temperature for 30 minutes.
- 9. Add 50 µl/well Stop Buffer (50 mM EDTA, pH 8) to stop the reaction.
- 10. Transfer 25  $\mu$ l of each reaction and 75  $\mu$ l dH<sub>2</sub>O/well to a 96-well streptavidincoated plate and incubate at room temperature for 60 minutes.
- 11. \*Wash three times with 200  $\mu\text{I/well PBS/T}$
- Dilute primary antibody, Phospho-Tyrosine mAb (P-Tyr-100), 1:1000 in PBS/T with 1% BSA. Add 100 μl/well primary antibody.
- **13.** Incubate at room temperature for 60 minutes.
- 14. \*Wash three times with 200 µl/well PBS/T
- For DELFIA<sup>®</sup> or Colorimetric ELISA detection methods please use the following protocols.

# **DELFIA®** Assay

- Prepare appropriate dilution of Europium labeled secondary antibody in PBS/T with 1% BSA (1:500 dilution for anti-mouse IgG or 1:1000 for anti-rabbit IgG).
- 2. Add 100  $\mu\text{I/well}$  secondary antibody solution.
- **3.** Incubate at room temperature for 30 minutes.
- 4. \*Wash five times with 200  $\mu\text{I/well PBS/T.}$
- 5. Add 100  $\mu\text{I/well}$  DELFIA® Enhancement Solution.
- 6. Incubate at room temperature for 5 minutes.
- Read plate using a Time Resolved Fluorescent plate reader using the following settings;
  - a. Excitation Filter: 340 nm
  - b. Emission Filter: 615 nm
  - c. Delay\*\*: 400 µs
- <sup>++</sup> Delay time is the delay from the excitation pulse to the beginning of the measurement.

## **Companion Products for DELFIA®**

DELFIA<sup>®</sup> Europium-labeled Anti-mouse IgG (PerkinElmer Life Sciences #AD0124) DELFIA<sup>®</sup> Europium-labeled Anti-rabbit IgG (PerkinElmer Life Sciences #AD0105) DELFIA<sup>®</sup> Enhancement Solution (PerkinElmer Life Sciences #1244-105) DELFIA<sup>®</sup> Streptavidin coated, 96-well, yellow plate (PerkinElmer Life Sciences AAAND-0005)

# **Colorimetric ELISA Assay**

- Prepare appropriate dilution of HRP labeled secondary antibody in PBS/T with 1% BSA (1:500 dilution for anti-mouse IgG or 1:1000 for anti-rabbit IgG).
- 2. Add 100 µl/well secondary antibody solution.
- 3. Incubate at room temperature for 30 minutes.
- **4.** \*Wash five times with 200 µl/well PBS/T.
- 5. Add 100 µl/well TMB substrate.
- 6. Incubate at room temperature for 15 minutes.
- 7. Add 100 µl/well of stop solution.
- 8. Mix well.
- 9. Read the absorbance at 450 nm with a microtiter plate reader.

#### **Companion Products For Colorimetric ELISA Assay**

Anti-mouse IgG, HRP Linked Antibody #7076 Anti-rabbit IgG, HRP Linked Antibody #7074 TMB Solution #7004 Stop Solution #7002

\* NOTE: Use of an automated microplate washer as well as centrifugation of plates when appropriate, greatly improves reproducibility.

Please contact Cell Signaling Technology for HTS-ready antibodies (PBS formulated and carrier-free), and detailed peptide substrate sequence information. Email: drugdiscovery@cellsignal.com