

#9656 Store at -20°C

StemLight™ Pluripotency Antibody Kit




1 Kit
 (100 tests)

Orders ■ 877-616-CELL (2355)
 orders@cellsignal.com
Support ■ 877-678-TECH (8324)
 info@cellsignal.com
Web ■ www.cellsignal.com

rev. 06/14/18

For Research Use Only. Not For Use In Diagnostic Procedures.

Products Included	Product #	Quantity	Application	Dilution	Isotype
Oct-4A (C30A3) Rabbit mAb	94310	100 tests	IF-IC	1:200	Rabbit IgG
Sox2 (D6D9) XP® Rabbit mAb	34516	100 tests	IF-IC	1:200	Rabbit IgG
Nanog (D73G4) XP® Rabbit mAb	99399	100 tests	IF-IC	1:200	Rabbit IgG
*SSEA4 (MC813) Mouse mAb	43782	100 tests	IF-IC	1:200	Mouse IgG3
*TRA-1-60(S) (TRA-1-60(S)) Mouse mAb	61220	100 tests	IF-IC	1:200	Mouse IgM
*TRA-1-81 (TRA-1-81) Mouse mAb	83321	100 tests	IF-IC	1:200	Mouse IgM



Projected confocal z-stack of human iPS cells using TRA-1-60(S) (TRA-1-60(S)) Mouse mAb (green, upper left), TRA-1-81 (TRA-1-81) Mouse mAb (green, upper middle), SSEA4 (MC813) Mouse mAb (green, upper right), Oct-4A (C30A3) Rabbit mAb (green, lower left), Sox2 (D6D9) XP® Rabbit mAb (green, lower middle) and Nanog (D73G4) XP® Rabbit mAb (green, lower right). Actin filaments were labeled with DyLight™ 554 Phalloidin #13054 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye).

Description: The StemLight™ Pluripotency Antibody Kit contains a panel of antibodies for the detection of proteins that are specifically expressed in human pluripotent cells. The kit can be used to track the pluripotent potential of human embryonic stem (ES) or induced pluripotent (iPS) cells. The loss of these markers indicates a loss of pluripotency or differentiation of the culture. The kit components are pre-optimized for parallel use in immunofluorescent analysis. Enough reagents are provided for 100 assays based on a working volume of 100 µl.

Specificity/Sensitivity: Each antibody in the StemLight™ Pluripotency Antibody Kit detects endogenous levels of their respective human pluripotency marker proteins.

Source/Purification: Nanog Antibody was produced by immunizing animals with a synthetic peptide corresponding to amino acid sequence at the amino terminus of human nanog. Antibodies are purified by Protein A and peptide affinity chromatography. Oct-4A antibody is produced by immunizing animals with recombinant protein specific to the amino terminus of human Oct-4A. Sox2 antibody is produced by immunizing animals with a synthetic peptide corresponding to amino acid sequences at the amino terminus of human Sox2. SSEA4, TRA-1-81, and TRA-1-60(S) antibodies are produced by immunizing animals with human embryonal carcinoma 2102Ep cl.2A6 cells.

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA and 50% glycerol. Store at -20°C. Do not aliquot the antibody.

**Note: Due to the extracellular location of the epitope, permeabilization (detergents or alcohols) should be omitted.*

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.

Background: Pluripotency is the ability of a cell to differentiate into cell types of the three germ layers, the endoderm, ectoderm and mesoderm. It is a property shared by embryonic stem cells, embryonic carcinoma and induced pluripotent cells.


Oct-4, Sox2 and Nanog are key transcriptional regulators that are highly expressed in pluripotent cells (1). Together they form a transcriptional network that maintains cells in a pluripotent state (2,3). Over-expression of Oct-4 and Sox2 along with Klf4 and c-Myc can induce pluripotency in both mouse and human somatic cells, highlighting their roles as key regulators of the transcriptional network necessary for renewal and pluripotency (4-5). It has also been demonstrated that overexpression of Oct-4, Sox2, Nanog and Lin28 can induce pluripotency in human somatic cells (6). Upon differentiation of pluripotent cultures, expression of Oct-4, Nanog and Sox2 is downregulated.

SSEA4, TRA-1-81 and TRA-1-60 antibodies recognize antigens expressed on the cell surface of all pluripotent cells. SSEA4 recognizes a glycolipid carbohydrate epitope (7). TRA-1-60(S) and TRA-1-81 antibodies recognize different proteoglycan epitopes on variants of the same protein, podocalyxin (8). These epitopes are neuraminidase sensitive and resistant, respectively. Reactivity of SSEA4, TRA-1-81 and TRA-1-60 antibodies with their respective cell surface markers are lost upon differentiation of pluripotent cells, corresponding with a loss of pluripotent potential (9).


U.S. Patent No. 5,675,063.
 DRAQ5 is a registered trademark of Biostatus Limited.
 DyLight is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

© 2015 Cell Signaling Technology, Inc. XP, StemLight and Cell Signaling Technology are trademarks of Cell Signaling Technology, Inc.


Applications Key: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry CHIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide
Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebra fish B—bovine
 Dg—dog Pg—pig Sc—S. cerevisiae All—all species expected Species enclosed in parentheses are predicted to react based on 100% sequence homology.



Confocal immunofluorescent analysis of NTERA-2 cells, untreated (left panel) or retinoic acid-treated (10 μ M all-trans RA for 5 days) (right panel), using TRA-1-60(S) (TRA-1-60(S)) Mouse mAb (green, upper), TRA-1-81 (TRA-1-81) Mouse mAb (green, middle) and SSEA4 (MC813) Mouse mAb (green, lower). Actin filaments have been labeled with DyLight™ 554 Phalloidin #13054 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye). Note the loss of pluripotency markers (green) as cells differentiate along the neuronal lineage with retinoic acid treatment.



Confocal immunofluorescent analysis of NTERA-2 cells, untreated (left panel) or retinoic acid-treated (10 μ M all-trans RA for 5 days) (right panel), using Neurofilament-L (C28E10) Rabbit mAb #2837 (green, upper), and a β -III tubulin antibody (green, lower). Actin filaments have been labeled with DyLight™ 554 Phalloidin #13054 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye). Note the appearance of neuronal markers and structures as cells differentiate along the neuronal lineage with retinoic acid treatment.



Confocal immunofluorescent analysis of NTERA-2 (left) and HeLa cells (right) using Oct-4A (C30A3) Rabbit mAb (green, upper), Sox2 (D6D9) XP® Rabbit mAb (green, middle) and Nanog (D73G4) XP® Rabbit mAb (green, lower).

Background References:

- (1) Looijenga, L.H. et al. (2003) *Cancer Res.* 63, 2244–2250.
- (2) Pesce, M. and Schöler, H.R. (2001) *Stem Cells* 19, 271–278.
- (3) Pan, G. and Thomson, J.A. (2007) *Cell. Res.* 17, 42–49.
- (4) Takahashi, K. and Yamanaka, S. (2006) *Cell* 126, 663–676.
- (5) Okita, K. et al. (2007) *Nature* 448, 313–317.
- (6) Yu, J. et al. (2007) *Science* 318, 1917–1920.
- (7) Henderson, J.K. et al. (2002) *Stem Cells* 20, 329–337.
- (8) Draper, J.S. et al. (2002) *J. Anat.* 200, 249–258.
- (9) Schopperle, W.M. and DeWolf, W.C. (2007) *Stem Cells* 25, 723–730.

Immunofluorescence Protocol

***IMPORTANT:** Please refer to the **APPLICATIONS** section on the front page of the data sheet to determine **IF THIS PRODUCT** is validated and approved for the specific protocol you will be using.

A Solutions and Reagents

NOTE: Prepare solutions with Milli-Q or equivalently purified water.

- 1. 10X Phosphate Buffered Saline (PBS):** To prepare 1 L add 80 g sodium chloride (NaCl), 2 g potassium chloride (KCl), 14.4 g sodium phosphate, dibasic (Na_2HPO_4) and 2.4 g potassium phosphate, monobasic (KH_2PO_4) to 1 L dH_2O . Adjust pH to 7.4.
- Formaldehyde, 16%, methanol free, Polysciences, Inc. (cat# 18814), use fresh, store opened vials at 4°C in dark, dilute in PBS for use.
- Xylene
- Ethanol, anhydrous denatured, histological grade, 100% and 95%
- Distilled water (dH_2O)
- 1X PBS/0.3% Triton X-100 (PBS/Triton):** To prepare 1 L, add 100 ml 10X PBS to 900 ml dH_2O . Add 3 ml Triton X-100 and mix.
- 10 mM Sodium Citrate Buffer:** To prepare 1 L, add 2.94 g sodium citrate trisodium salt dihydrate ($\text{C}_6\text{H}_5\text{Na}_3\text{O}_7 \cdot 2\text{H}_2\text{O}$) to 1 L dH_2O . Adjust pH to 6.0.
- 1X PBS, high salt (0.4M) (high salt PBS):** To prepare 1L, add 100 ml 10X PBS to 900 ml dH_2O . Add 23.38 g NaCl and mix.
- Fluorochrome-conjugated secondary antibody

NOTE: When using any primary or fluorochrome-conjugated secondary antibody for the first time, titrate the antibody to determine which dilution allows for the strongest specific signal with the least background for your sample.

10. Prolong[®] Gold Antifade Reagent (Invitrogen, Eugene, OR, Cat# P36930)

B Specimen Preparation

I. Cultured Cell Lines (IF-IC)

IMPORTANT: Please check the **APPLICATIONS** section of the data sheet to verify that this product is validated and approved for **(IF-IC)**.

NOTE: This general fixation protocol will work with most antibodies and cell lines. However, we recommend you try different IF/IC fixation methods (methanol or acetone alone, aldehyde alone, or combinations of these) to identify the optimal fixation protocol for each antibody and/or cell line.

NOTE: Cells should be grown, treated, fixed, and stained directly in multiwell plates, chamber slides, or on coverslips.

1. Rinse cells briefly in PBS.
2. Aspirate PBS, cover cells to a depth of 2-3 mm with 2-4% formaldehyde in PBS.

NOTE: Formaldehyde is toxic, use only in fume hood.

3. Allow cells to fix for 15 minutes at room temperature.
4. Aspirate fixative, rinse three times in PBS for 5 minutes each.

5. Methanol Permeabilization Step (if required, please refer to front page): After formaldehyde fixation, cover cells with ice-cold 100% methanol (use enough to cover cells completely to a depth of 3-5 mm, **DO NOT LET CELLS DRY**), incubate cells in methanol for 10 minutes in freezer, rinse in PBS for 5 minutes.

6. Proceed with Immunostaining section C.

II. Paraffin Sections (IF-P)

IMPORTANT: Please check the **APPLICATIONS** section of the data sheet to verify that this product is validated and approved for **(IF-P)**.

Deparaffinization/Rehydration:

1. Incubate sections in three washes of xylene for 5 minutes each.
2. Incubate sections in two washes of 100% ethanol for 10 minutes each.
3. Incubate sections in two washes of 95% ethanol for 10 minutes each.
4. Rinse sections twice in dH_2O for 5 minutes each.

Antigen Unmasking:

1. Place slides in room temperature 10 mM sodium citrate buffer pH 6.0.
2. Bring slides to boiling in sodium citrate buffer using water bath or microwave, then maintain at 95-99°C for 10 minutes.
3. Cool slides for 30 minutes on bench top.
4. Rinse sections in dH_2O three times for 5 minutes each.
5. Rinse sections in PBS for 5 minutes.
6. Proceed with Immunostaining section C.

III. Frozen/Cryostat Sections (IF-F)

IMPORTANT: Please check the **APPLICATIONS** section of the data sheet to verify that this product is validated and approved for **(IF-F)**.

NOTE: Fresh frozen/unfixed sections should be fixed immediately in 2-4% formaldehyde as follows to preserve signaling epitopes.

1. Cover sections with 2-4% formaldehyde in PBS

NOTE: Formaldehyde is toxic, use only in fume hood.

2. Allow cells to fix for 15 minutes at room temperature.
3. Rinse slides three times in PBS for 5 minutes each.

C Immunostaining

NOTE: All subsequent incubations should be carried out at room temperature unless otherwise noted in a humid light-tight box or covered dish/plate to prevent drying and fluorochrome fading.

1. Block specimen in 5% normal serum from same species as secondary antibody (eg. normal goat serum, normal donkey serum) in PBS/Triton for 60 minutes.
2. While blocking, prepare primary antibody by diluting as indicated on datasheet in PBS/Triton. You will need 50-100 μl per section, 25-50 μl per coverslip, chamber, or well (48 or 96 well plate).
3. Aspirate blocking solution, apply diluted primary antibody.

NOTE: For double-labeling, prepare a cocktail of mouse and rabbit primary antibodies at their appropriate dilutions in PBS/Triton.

4. Incubate overnight at 4°C.
5. Rinse three times in PBS for 5 minutes each.

OPTION: To decrease background stain, rinse in high salt PBS for two minutes between second and third PBS rinses. Be aware, this may reduce specific staining of some antibodies.

NOTE: If using primary antibodies directly conjugated with AlexaFluor[®] fluorochromes, then skip to step C8.

6. Incubate in fluorochrome-conjugated secondary antibody diluted in PBS/Triton for 1-2 hours at room temperature in dark.

NOTE: For double-labeling, prepare a cocktail of fluorochrome-conjugated anti-mouse and anti-rabbit primary antibodies at their appropriate dilutions in PBS/Triton.

7. Rinse in PBS/high salt PBS as in step 5.
8. Coverslip slides with Prolong[®] Gold Antifade Reagent or apply just enough to cover cells in multiwell plate.
9. Seal slides by painting around edges of coverslips with nail polish.
10. Examine specimens immediately using appropriate excitation wavelength, depending on fluorochrome for best results or store flat at 4°C in dark.