

# Caspase-7 Activity Assay Kit

✓ 1 Kit (40 assays)

rev02/04



Cell Signaling  
TECHNOLOGY®

Orders ■ 877-616-CELL (2355)  
orders@cellsignal.com

Support ■ 877-678-TECH (8324)  
info@cellsignal.com

Web ■ www.cellsignal.com

**Description:** The Caspase-7 Activity Assay includes all reagents required to measure Caspase-7 activity, including antibody, substrate, reaction buffer and cell extracts buffer. The assay system uses Caspase-7 antibody to selectively immunoprecipitate Caspase-7. Fluorometric substrate is added and is cleaved proportionally to the amount of activated Caspase-7, generating free fluorescent AFC. Free AFCs are determined fluorometrically at  $\lambda_{max} = 505$  nm. Compared to other commercially available Caspase-7 activity assays, this kit is designed to be more specific by using Caspase-7 antibody to selectively immunoprecipitate Caspase-7 from cell lysates.

**Species Cross Reactivity:** H

#### Kit Components:

**Caspase-7 Antibody:** Polyclonal antibodies are produced by immunizing rabbits with a synthetic peptide (KLH coupled) corresponding to human Caspase-7. Antibodies are purified by protein A and peptide affinity chromatography and supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA and 50% glycerol. Store at -20°C.

#### Fluorogenic Caspase-7 Substrate

Ac-Asp-Glu-Val-Asp-AFC

**10X Cell Extracts Buffer:** 1X concentration: 50mM Pipes/KOH (pH 6.5), 2mM EDTA, 0.1% Chaps, 5mM DTT, 20µg/ml Leupeptin, 10µg/ml pepstatin A, 10 µg/ml aprotinin

#### 10X Reaction/Washing Buffer

1X concentration: 50mM HEPES (pH 7.5), 75mM NaCl, 5mM DTT, 0.1% Chaps 200X DTT (1 M DTT) 1 ml

**Background:** Caspase-7 (CMH-1, Mch3, ICE-LAP3) has been identified as a major contributor to the execution of apoptosis (2-4). Caspase-7 is an effector caspase (along with caspase-2 and -3), meaning that it cleaves essential cellular machinery rather than activating other caspases (5-8). Caspase-7 is cleaved by many enzymes, including caspases-3, -6, -8, -9 and granzyme B (1,4,5). Once activated, caspase-7 cleaves many of the same substrates as caspase-3, including poly (ADP-ribose) polymerase, or PARP (2,4).

#### Background References:

(1) Nuñez, G. et al. (1998) *Oncogene* 17, 3237-3245.

#### Products Included in Kit

Products	Product #	Kit Quantity
Caspase-7 Antibody	[Antibody is unique to this kit.]	150 µl
Fluorogenic Caspase-7 Substrate		100 µl
Cell Extract Buffer (10X)	9852	5 ml
10X Reaction Buffer		15 ml
200X Dithiothreitol (1 M DTT)		1 ml

## Caspase-7 Assay Kit Overview


### Step 1: Selective IP of Caspase-7 using Caspase-7 Antibody

- Add Caspase-7 Antibody
- Add protein A sepharose beads\*

### Step 2: Incubate IP pellets in Reaction Buffer containing Caspase-7 Substrate (Ac-DEVD-AFC). Active Caspase-7 cleaves Ac-DEVD-AFC between DEVD and AFC, generating free fluorescent AFC.

### Step 3: Measure free fluorescent AFC at excitation 400nm, emission 505nm.

\* reagent not included in kit




Titration curve of Caspase-7 activity assay using varying amounts of Jurkat cells. Cells were either treated with etoposide for 6 hours or left untreated. Cell extracts were prepared for immunoprecipitation using Caspase-7 Antibody. Immunocomplexes were incubated with peptide substrate Ac-DEVD-AFC. Specific Caspase-7 activity was measured based on the fluorescent AFC group from the hydrolysis of Ac-DEVD-AFC by activated Caspase-7. Specific inhibition assay was performed by preincubating the immunocomplex with Ac-DEVD-CHO for 30 minutes before adding Ac-DEVD-AFC.

- (2) Fernandes-Alnemri, T. et al. (1995) *Cancer Res.* 55, 6045–6052.
- (3) Duan, H. et al. (1996) *J. Biol. Chem.* 271, 1621–1625.
- (4) Lippke, J.A. et al. (1996) *J. Biol. Chem.* 271, 1825–1828.
- (5) Cohen, G.M. (1997) *Biochem. J.* 326, 1–16.
- (6) Thornberry, N.A. et al. (1997) *J. Biol. Chem.* 272, 17907–17911.
- (7) Chandler, J.M. et al. (1998) *J. Biol. Chem.* 273, 10815–10818.
- (8) MacFarlane, M. et al. (1997) *J. Cell Biol.* 137, 469–479.

**Companion Products:**

- Cleaved Caspase-7 (Asp198) Antibody #9491
- Caspase-7 Antibody #9492
- Chaps Cell Exrtact Buffer (10x) #9852
- Caspase-3 (3G2) mAb (human specific) #9668
- Cleaved Caspase-3 (Asp175) (5A1) Rabbit Monoclonal Antibody #9664
- Caspase-3 (10G1) Rabbit Monoclonal Antibody #9665



Recommended Antibody Dilutions	
	Caspase-7 Antibody
Immunoprecipitation:	1:100

Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken X—Xenopus Z—zebra fish All—all species expected  
Species enclosed in parentheses are predicted to react based on 100% sequence homology.

Orders ■ 877-616-CELL (2355) orders@cellsignal.com Support ■ 877-678-TECH (8324) info@cellsignal.com Web ■ www.cellsignal.com

## Immunoprecipitation/Cleavage Reaction Protocols

### A Solutions and Reagents

**Note:** Prepare solutions with Milli-Q or equivalently purified water.

- A1 Cell Extracts Buffer:** Immediately before use add DTT to final concentration of 5mM (1:200 dilution) and PMSF\* to final concentration to 1mM.
- A2 Reaction/Washing Buffer:** immediately before use add DTT to final concentration of 5mM (1:200 dilution).
- A3 \*Protein A Sepharose Beads:** (can be stored for 2 weeks at 4°C) Add 5 ml of 1X PBS to 1.5g of protein A sepharose beads, shake for 2 hours at 4°C. Spin down, wash pellet twice with PBS and resuspend beads in 1 volume of PBS.
- \*reagent not included in kit

### B Preparing Cell Lysates

- B1** Treat cells with regulator for desired time.
- B2** To harvest cells under **nondenaturing conditions**, remove media and rinse cells once with ice-cold PBS.
- B3** Pellet the cells and resuspend them in the appropriate volume of 1X ice-cold Cell Extracts Buffer plus 1mM PMSF.
- B4** Freeze/thaw cells 3 times.
- B5** Microcentrifuge at top speed for 10 minutes at 4°C; transfer the supernatant to a new tube. The supernatant is the cell lysate. If necessary, lysate can be stored at -80°C.

### C Immunoprecipitation using Caspase 7 Antibodies:

- C1** To 300  $\mu$ l cell lysate, add Caspase-7 Antibody (1:100 dilution); incubate with gentle rocking for 2 hours.
- C2** Add protein A sepharose beads (10-20  $\mu$ l [50% beads]). Incubate with gentle rocking for 1-3 hours at 4°C.
- C3** Microcentrifuge at 4,000 rpm for 30 seconds at 4°C. Remove unbound supernatant. Wash pellet twice with 1 ml of 1X Reaction Buffer. Keep on ice.

### D Cleavage Reaction and Detection:

- D1** Add 100  $\mu$ L of Reaction Buffer containing a final concentration of 0.4mM Ac-DEVD-AFC (2  $\mu$ L of the substrate stock solution) to the bead pellet.
- D2** Incubate at 37°C overnight.
- D3** The supernatant contains the free AFC. Measure AFC at excitation 400nm, emission 505nm.