

#7002 Store at +4°C

STOP Solution

P4: 50 ml (2 x 25 ml)
P6: 250 ml (10 x 25 ml)



Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Description

STOP Solution is a proprietary, ready to use from the bottle solution, used to terminate the peroxidase/TMB reaction for ELISA applications. The TMB substrate reacts with immobilized horseradish peroxidase (HRP) conjugated secondary antibodies to produce a blue solution. Color intensity is an indication of analyte level. After attaining the desired intensity, the reaction is terminated by addition of STOP Solution. Upon addition of STOP Solution the color turns from blue to yellow. Absorbance at 450 nm can be read immediately, and the reaction is stable for one hour.

Storage

Store at 4°C.

Background

STOP Solution is used to terminate TMB/peroxidase reaction for ELISA applications and is intended for use in combination with CST's PathScan® Sandwich ELISA Antibody Pairs.

Directions for Use

- (a) After incubation with horseradish peroxidase (HRP) conjugated secondary antibody, wash plates with PBS/0.05% Tween-20.
- (b) Add 100 µl of TMB solution/well. Incubate for a specific amount of time or until desired signal is achieved. See product datasheet for recommended incubation times for specific antibody pair.
- (c) Add 100 µl of STOP Solution.
- (d) Read at Absorbance 450 nm.

Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.