

Pazopanib



Orders: 877-616-CELL (2355) orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com

cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Background

Pazopanib is a multikinase inhibitor that potently targets VEGFR1 (IC $_{50}$ = 10 nM), VEGFR2 (IC $_{50}$ = 30 nM), VEGFR3 (IC $_{50}$ = 47 nM), PDGFR α (IC $_{50}$ = 71 nM), PDGFR β (IC $_{50}$ = 84 nM), and c-Kit (IC $_{50}$ = 74 nM) tyrosine kinases involved in tumor progression and angiogenesis, and can also inhibit many other tyrosine kinases at nanomolar concentrations (1). Research studies have demonstrated that pazopanib effectively blocks ligand-induced autophosphorylation of VEGFR2, PDGFR β , and c-Kit *in vitro* (1,2), and selectively inhibits VEGF-induced HUVE cell proliferation over FGF (IC $_{50}$ = ~21 nM vs ~720 nM). Investigators have demonstrated that pazopanib inhibits the growth, survival, and migration of multiple myeloma (MM) cell types (3).

Molecular Formula C₂₁H₂₃N₇O₂S Molecular Weight 437.52 g/mol

Purity >99%

CAS 444731-52-6

Solubility Soluble in DMSO at 8.3mg/ml.

Solubliney Soluble in Bings at olding in.

Store lyophilized or in solution at -20°C, desiccated. In lyophilized form, the chemical is stable for 24 months. Once in solution, use within 3 months to prevent loss of potency. Aliquot to avoid multiple

freeze/thaw cycles.

Directions for Use

Storage

Pazopanib is supplied as a lyophilized powder. For a 10 mM stock, reconstitute the 10 mg in 2.29 ml DMSO. Working concentrations and length of treatment can vary depending on the desired effect, but it is typically used as a pretreatment at 0.1-10 μ M for 0.5-2 hr prior to treating with a stimulator. It can also be used alone, with varying treatment times lasting up to 24 hr.

Solubility: Soluble in DMSO at 8 mg/mL with slight warming; very poorly soluble in ethanol and water with maximum in water \sim 10-20 μ M.

Background References

- 1. Kumar, R. et al. (2007) Mol Cancer Ther 6, 2012-21.
- 2. Kumar, R. et al. (2009) Br J Cancer 101, 1717-23.
- 3. Podar, K. et al. (2006) Proc Natl Acad Sci U S A 103, 19478-83.

Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

XP is a registered trademark of Cell Signaling Technology, Inc.

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any

commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.