Revision 2

SM-164	Cell Signaling
	Orders: 877-616-CELL (2355) orders@cellsignal.com
2 1 mg	Support: 877-678-TECH (8324)
1 mg 60095#	Web: info@cellsignal.com cellsignal.com
#2	3 Trask Lane Danvers Massachusetts 01923 USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Background	SM-164 is a bivalent, cell-permeable, non-peptide, SMAC mimetic and potent XIAP inhibitor with an IC ₅₀ value of 1.39 nM. SM-164 binds to the BIR2 and BIR3 domains of XIAP and effectively inhibits cell growth and induces apoptosis in HL-60 leukemia cells (1). It also activates caspase-3 and caspase-9 that, in turn, induce apoptosis in leukemia cells (1,2). In breast cancer mouse xenograft models, SM-164 in combination with TNF-related apoptosis-inducing ligand (TRAIL) has been shown to reduce tumor volume by 80%. SM-164 is an effective apoptotic inducer of tumor cells and tissues and is a promising candidate for a new class of anticancer drugs (2). Treatment of RIP3-positive cells with SM-164 in combination with TNF-alpha and the caspase inhibitor Z-VAD, leads to necroptosis (3).
Molecular Formula	C ₆₂ H ₈₄ N ₁₄ O ₆
Molecular Weight	1121.4 g/mol
Purity	>98%
CAS	957135-43-2
Solubility	Soluble in DMSO at 50 mg/mL.
Storage	Store lyophilized at -20°C, desiccated. In lyophilized form, the chemical is stable for 24 months. Once in solution, store at -20°C and use within 3 months to prevent loss of potency. <i>Aliquot to avoid multiple freeze/thaw cycles.</i>
Directions for Use	SM-164 is supplied as a lyophilized powder. For a 10 mM stock, reconstitute 1 mg of powder in 89 μ L of DMSO. Working concentrations and length of treatment can vary depending on the desired effect.
Background References	1. Sun, H. et al. (2007) <i>J Am Chem Soc</i> 129, 15279-94. 2. Lu, J. et al. (2008) <i>Cancer Res</i> 68, 9384-93. 3. Sun, L. et al. (2012) <i>Cell</i> 148, 213-27.
Trademarks and Patents	Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.
	All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.
	Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.
	Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with

Orders: 877-616-CELL (2355) • orders@cellsignal.com • Support: 877-678-TECH (8324) • info@cellsignal.com • Web: cellsignal.com For Research Use Only. Not for Use in Diagnostic Procedures.