## Tris Buffered Saline with Tween® 20 (TBST - 10X)



Supplied as a 10X stock. **1000 ml** 



Orders 877-616-CELL (2355) orders@cellsignal.com Support 877-678-TECH (8324) info@cellsignal.com Web www.cellsignal.com

rev. 08/22/18

## For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Tris buffered saline (TBS) solution with the detergent Tween® 20 for use as a wash buffer during western blotting. Product is shipped and stored at room temperature. 1X Formulation: 137 mM Sodium Chloride, 20 mM Tris, 0.1% Tween® 20. Supplied at pH 7.6.

Directions for Use: Dilute Tris Buffered Saline with Tween® 20 (TBST-10X) to a 1X solution using ddH<sub>2</sub>O. This product supplies enough 10X material to make 10 liters of 1X solution.

Storage: Supplied as a 10X solution. Store at room temperature.

Please visit www.cellsignal.com for a complete listing of recommended companion products.