

#1325

Store at -20°C

PTP1B (Tyr66) Biotinylated peptide

✓ 1.25 ml at 6µM



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 03/21/16

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: This biotinylated peptide contains the residues surrounding Tyr66 of PTP1B. It was generated for use in CST's HTScan® kinase assay kits, but may also serve as a substrate in other heterogeneous or homogeneous kinase assays.

Peptide Core Sequence: DNDY*INA

Molecular Weight: 2141 daltons

Quality Control: The quality of the biotinylated peptide was evaluated by reverse-phase HPLC and by mass spectrometry.

Directions for Use: The phosphorylated form of the peptide can be detected with the Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411. Sample kinase assay protocols can be found on corresponding kinase assay kit data sheets (see Companion Products).

Storage: Supplied in 0.0001% DMSO. Store at -20°C.

Companion Products:

HTScan® EGF Receptor Kinase Assay Kit #7909

Tyrosine Kinase Substrate Screening Kit #7450

Phospho-Tyrosine Mouse mAb (P-Tyr-100) #9411