

#14282 Store at 4°C

SimpleChIP® Enzymatic Cell Lysis Buffers A & B

www.cellsignal.com

1 Kit (30 immunoprecipitations)

Support: 877-678-TECH (8324)
www.cellsignal.com/support

Orders: 877-616-CELL (2355)
orders@cellsignal.com

rev. 04/10/17

For Research Use Only. Not For Use In Diagnostic Procedures.

Products Included	Product #	Quantity
Buffer A (4X)	7006	25 ml
Buffer B (4X)	7007	25 ml

Description: This product is offered to conveniently provide additional Buffer A and Buffer B reagents for harvesting cells and tissues and preparing chromatin using our SimpleChIP® (#9002, #9003) and SimpleChIP® Plus (#9004, #9005) Enzymatic Chromatin IP Kits. These kits provide all the reagents required for performing 6 chromatin preparations (or optimizations) and 30 chromatin immunoprecipitation (ChIP) assays, however there are instances where extra Buffer A and Buffer B reagents are desired.

Storage: Store both buffers at 4°C. They are stable for 12 months upon receipt when stored properly.

Directions for Use: Use as cell or tissue harvesting and chromatin preparation buffers in chromatin immunoprecipitation (ChIP) assays as directed in the protocols for the SimpleChIP® and SimpleChIP® Plus Enzymatic Chromatin IP Kits.

For product specific protocols and a complete listing of recommended companion products please see the product web page at www.cellsignal.com

Thank you for your recent purchase. If you would like to provide a review visit www.cellsignal.com/comments.

© 2014 Cell Signaling Technology, Inc.
SimpleChIP and Cell Signaling Technology are trademarks of Cell Signaling Technology, Inc.



Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide **Species Cross-Reactivity:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected *Species enclosed in parentheses are predicted to react based on 100% homology.