



Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

Store at +4C
#14791

Phospho-Zap-70 (Tyr319)/Syk (Tyr352) (65E4) Rabbit mAb (PE Conjugate)

For Research Use Only. Not for Use in Diagnostic Procedures.

Applications: FC-FP	Reactivity: H M	Sensitivity: Endogenous	Source/Isotype: Rabbit IgG	UniProt ID: #P43403, #P43405	Entrez-Gene Id: 7535, 6850
-------------------------------	---------------------------	-----------------------------------	--------------------------------------	--	--------------------------------------

Product Usage Information	Application	Dilution
Storage	Flow Cytometry (Fixed/Permeabilized)	1:50
Specificity/Sensitivity	Supplied in PBS (pH 7.2), less than 0.1% sodium azide and 2 mg/ml BSA. Store at 4°C. Do not aliquot the antibodies. Protect from light. Do not freeze.	
Species predicted to react based on 100% sequence homology	Phospho-Zap-70 (Tyr319)/Syk (Tyr352) (65E4) Rabbit mAb (PE Conjugate) detects endogenous levels of Zap-70 only when phosphorylated at Tyr319. It cross-reacts with endogenous levels of Syk when phosphorylated at Tyr352.	
Source / Purification	Rat, Hamster, Chicken, Bovine, Dog, Pig, Horse	
Description	Monoclonal antibody is produced by immunizing animals with a synthetic phosphopeptide corresponding to residues surrounding Tyr319 of human Zap-70 protein.	
Background	This Cell Signaling Technology antibody is conjugated to phycoerythrin (PE) and tested in-house for direct flow cytometry analysis in human cells. The antibody is expected to exhibit the same species cross-reactivity as the unconjugated Phospho-Zap-70 (Tyr319)/Syk (Tyr352) (65E4) Rabbit mAb #2717.	
Background References	The Syk family protein tyrosine kinase Zap-70 is expressed in T and NK cells and plays a critical role in mediating T cell activation in response to T cell receptor (TCR) engagement (1). Following TCR engagement, Zap-70 is rapidly phosphorylated on several tyrosine residues through autophosphorylation and transphosphorylation by the Src family tyrosine kinase Lck (2-6). Tyrosine phosphorylation correlates with increased Zap-70 kinase activity and downstream signaling events. Expression of Zap-70 is correlated with disease progression and survival in patients with chronic lymphocytic leukemia (7,8).	
	<ol style="list-style-type: none"> 1. Chu, D.H. et al. (1998) <i>Immunol Rev</i> 165, 167-80. 2. Iwashima, M. et al. (1994) <i>Science</i> 263, 1136-9. 3. Neumeister, E.N. et al. (1995) <i>Mol Cell Biol</i> 15, 3171-8. 4. Chan, A.C. et al. (1995) <i>EMBO J</i> 14, 2499-508. 5. Williams, B.L. et al. (1999) <i>EMBO J</i> 18, 1832-44. 6. Di Bartolo, V. et al. (1999) <i>J Biol Chem</i> 274, 6285-94. 7. Wiestner, A. et al. (2003) <i>Blood</i> 101, 4944-51. 8. Crespo, M. et al. (2003) <i>N Engl J Med</i> 348, 1764-75. 	

Species Reactivity Species reactivity is determined by testing in at least one approved application (e.g., western blot).

Applications Key **FC-FP:** Flow Cytometry (Fixed/Permeabilized)

Cross-Reactivity Key **H:** Human **M:** Mouse

Trademarks and Patents Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.
All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no

force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.