



Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

Store at -20C
#4896

PLK3 (D14F12) Rabbit mAb

For Research Use Only. Not for Use in Diagnostic Procedures.

Applications: W	Reactivity: H M	Sensitivity: Endogenous	MW (kDa): 80	Source/Isotype: Rabbit IgG	UniProt ID: #Q9H4B4	Entrez-Gene Id: 1263
---------------------------	---------------------------	-----------------------------------	------------------------	--------------------------------------	-------------------------------	--------------------------------

Product Usage Information

Application

Western Blotting

Dilution

1:1000

Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

Specificity/Sensitivity

PLK3 (D14F12) Rabbit mAb detects endogenous levels of total PLK3 protein.

Species predicted to react based on 100% sequence homology

Rat, Monkey

Source / Purification

Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Cys625 of human PLK3 protein.

Background

At least 4 distinct polo-like kinases exist in mammalian cells: PLK1, PLK2, PLK3 and PLK4/SAK (1). Like the other PLK family members, PLK3 contains an amino-terminal catalytic domain and a conserved carboxy-terminal domain termed the Polo box (2). PLK3, also called proliferation-related kinase (Prk) (3), was originally described as a fibroblast growth factor (FGF)-inducible kinase (Fnk) and identified as an immediate-early response gene responsive to FGF-1 and other mitogens (4). PLK3 is a cytokine-inducible serine/threonine kinase whose protein expression is cell cycle regulated. Though its expression is found primarily in G1 phase of the cell cycle, PLK3 is detected in G0 and in late telophase prior to cytokinesis (5). Like the other PLK family members, PLK3 functions mainly as a regulator of the cell cycle. Specifically, PLK3 is required for entry into S phase and is a critical regulator of G1 events, as indicated by RNAi-induced PLK3-depleted cells (2). PLK3 is also involved in the regulation of DNA damage response via phosphorylation of p53 on Ser20 (6). PLK3 may act as a tumor suppressor as Plk3-deficient mice develop spontaneous tumors in various organs (7). Unlike PLK1, PLK3 expression is down regulated in cancers including lung (3), head and neck (8), and colon (9).

Background References

1. Nigg, E.A. (1998) <cite>Curr Opin Cell Biol</cite> 10, 776-83.
2. Zimmerman, W.C. and Erikson, R.L. (2007) <cite>Proc Natl Acad Sci USA</cite> 104, 1847-52.
3. Li, B. et al. (1996) <cite>J Biol Chem</cite> 271, 19402-8.
4. Donohue, P.J. et al. (1995) <cite>J Biol Chem</cite> 270, 10351-7.
5. Zimmerman, W.C. and Erikson, R.L. (2007) <cite>Cell Cycle</cite> 6, 1314-8.
6. Xie, S. et al. (2001) <cite>J Biol Chem</cite> 276, 43305-12.
7. Yang, Y. et al. (2008) <cite>Cancer Res</cite> 68, 4077-85.
8. Dai, W. et al. (2000) <cite>Genes Chromosomes Cancer</cite> 27, 332-6.
9. Dai, W. et al. (2002) <cite>Int J Oncol</cite> 20, 121-6.

Species Reactivity

Species reactivity is determined by testing in at least one approved application (e.g., western blot).

Western Blot Buffer

IMPORTANT: For western blots, incubate membrane with diluted primary antibody in 5% w/v BSA, 1X TBS, 0.1% Tween® 20 at 4°C with gentle shaking, overnight.

Applications Key

W: Western Blotting

Cross-Reactivity Key

H: Human **M:** Mouse

Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.