

MLK1 (D66E2) Rabbit mAb

Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Applications: W, IP	Reactivity: H M	Sensitivity: Endogenous	MW (kDa): 120	Source/Isotype: Rabbit IgG	UniProt ID: #P80192	Entrez-Gene Id: 4293
-------------------------------	---------------------------	-----------------------------------	-------------------------	--------------------------------------	-------------------------------	--------------------------------

Product Usage Information	Application Western Blotting Immunoprecipitation	Dilution 1:1000 1:50
Storage	Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.	
Specificity/Sensitivity	MLK1 (D66E2) Rabbit mAb detects endogenous levels of total MLK1 protein. The antibody also cross-reacts with a 45 kDa protein of unknown origin.	
Source / Purification	Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Leu1070 of human MLK1 protein.	
Background	Mixed-lineage kinases (MLKs) belong to the mitogen activated kinase kinase kinase (MAPKKK) family of dual-specificity protein kinases. While not particularly well conserved at the sequence level, MLK1, 2 and 3 share a conserved domain structure consisting of a catalytic core and two isoleucine/leucine zipper motifs among other protein-protein binding domains (1). MLK1 preferentially stimulates the JNK (c-Jun amino-terminal kinase) pathway in response to agonists and stress (2). Although multiple phosphorylation events are required for full activation of MLK1, two autophosphorylation sites within the activation loop (Ser308 and Thr312) appear to be the predominant activation residues (3). In neuronal cells, MLK1 appears to function downstream of the small G-proteins Rac1 and Cdc42 and upstream of MKK4 and MKK7 to promote apoptosis (2).	
Background References	<ol style="list-style-type: none"> 1. Dorow, D.S. et al. (1993) <i>Eur J Biochem</i> 213, 701-10. 2. Xu, Z. et al. (2001) <i>Mol Cell Biol</i> 21, 4713-24. 3. Durkin, J.T. et al. (2004) <i>Biochemistry</i> 43, 16348-55. 	

Species Reactivity	Species reactivity is determined by testing in at least one approved application (e.g., western blot).
Western Blot Buffer	IMPORTANT: For western blots, incubate membrane with diluted primary antibody in 5% w/v BSA, 1X TBS, 0.1% Tween@ 20 at 4°C with gentle shaking, overnight.
Applications Key	W: Western Blotting IP: Immunoprecipitation
Cross-Reactivity Key	H: Human M: Mouse
Trademarks and Patents	Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc. All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.
Limited Uses	<p>Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.</p> <p>Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a</p>

component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.