

#5425

Store at -20°C

Normal Goat Serum

10 ml



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 11/29/18

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Premium quality normal goat serum for use as blocking reagent and antibody diluent for chromogenic and fluorescent immunohistochemical and immunocytochemical assays. Typically, normal serum from the same species as the secondary antibody is used in the blocking and diluent buffers.

Storage: This product contains less than 0.1% of thimerosal and is stable when stored at 4°C for up to one month. For longer term storage, aliquot into working volumes and store at -20°C or -80°C.

Please visit www.cellsignal.com for a complete listing of recommended companion products.