

#9804

Store at -20°C

ATP (10 mM)

1 ml
(100 Assays (96 Well Format))



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 02/03/17

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Adenosine-5'-triphosphate (ATP) supplied as a 10 mM solution in doubly distilled water as a disodium salt.

Directions for Use: For *in vitro* kinase assays, we recommend a final ATP concentration of 200 μ M.

Companion Products:

Kinase Buffer (10X) #9802

Cell Lysis Buffer (10X) #9803