

Akt (pan) Blocking Peptide

✓ 100 µg



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com

rev. 09/13/18

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: This peptide is used to block Akt (pan) (11E7) Rabbit mAb #4685 reactivity.

Background: Akt, also referred to as PKB or Rac, plays a critical role in controlling survival and apoptosis (1-3). This protein kinase is activated by insulin and various growth and survival factors to function in a wortmannin-sensitive pathway involving PI3 kinase (2,3). Akt is activated by phospholipid binding and activation loop phosphorylation at Thr308 by PDK1 (4) and by phosphorylation within the carboxy terminus at Ser473. The previously elusive PDK2 responsible for phosphorylation of Akt at Ser473 has been identified as mammalian target of rapamycin (mTOR) in a rapamycin-insensitive complex with rictor and Sin1 (5,6). Akt promotes cell survival by inhibiting apoptosis through phosphorylation and inactivation of several targets, including Bad (7), forkhead transcription factors (8), c-Raf (9), and caspase-9. PTEN phosphatase is a major negative regulator of the PI3 kinase/Akt signaling pathway (10). LY294002 is a specific PI3 kinase inhibitor (11). Another essential Akt function is the regulation of glycogen synthesis through phosphorylation and inactivation of GSK-3 α and β (12,13). Akt may also play a role in insulin stimulation of glucose transport (12). In addition to its role in survival and glycogen synthesis, Akt is involved in cell cycle regulation by preventing GSK-3 β -mediated phosphorylation and degradation of cyclin D1 (14) and by negatively regulating the cyclin dependent kinase inhibitors p27 Kip1 (15) and p21 Waf1/Cip1 (16). Akt also plays a critical role in cell growth by directly phosphorylating mTOR in a rapamycin-sensitive complex containing raptor (17). More importantly, Akt phosphorylates and inactivates tuberin (TSC2), an inhibitor of mTOR within the mTOR-raptor complex (18,19).

Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. The peptide detects Akt (pan) (11E7) Rabbit mAb #4685 by peptide dot blot.

Directions for Use: Use as a blocking reagent to evaluate the specificity of antibody reactivity in peptide dot blot protocols.

Background References:

- (1) Franke, T.F. (1997) *Cell* 88, 435–437.
- (2) Burgering, B.T. and Coffey, P.J. (1995) *Nature* 376, 599–602.
- (3) Franke, T.F. et al. (1995) *Cell* 81, 727–736.
- (4) Alessi, D.R. et al. (1996) *EMBO J.* 15, 6541–6551.
- (5) Sarbassov, D.D. et al. (2005) *Science* 307, 1098–1101.
- (6) Jacinto, E. et al. (2006) *Cell* 127, 125–137.
- (7) Cardone, M.H. et al. (1998) *Science* 282, 1318–1321.
- (8) Brunet, A. et al. (1999) *Cell* 96, 857–868.
- (9) Zimmerman, S. et al. (1999) *Science* 286, 1741–1744.
- (10) Cantley, L.C. et al. (1999) *Proc. Natl. Acad. Sci. USA* 96, 4240–4245.
- (11) Vlahos, C. et al. (1994) *J. Biol. Chem.* 269, 5241–5248.
- (12) Hajduch, E. et al. (2000) *FEBS Lett.* 492, 199–203.
- (13) Cross, D.A. et al. (1995) *Nature* 373, 785–789.
- (14) Diehl, J.A. et al. (1998) *Genes Dev.* 12, 3499–3511.
- (15) Gesbert, F. et al. (2000) *J. Biol. Chem.* 275, 39223–39230.
- (16) Zhou, B.P. et al. (2001) *Nat. Cell Biol.* 3, 245–252.
- (17) Nave, B.T. et al. (1999) *Biochem. J.* 344, 427–431.
- (18) Manning, B.D. et al. (2000) *Mol. Cell* 4, 648–657.
- (19) Manning, B.D. et al. (2002) *Mol. Cell* 10, 151–162.
- (20) Inoki, K. et al. (2002) *Nat. Cell Biol.* 4, 648–657.

Entrez-Gene ID #207

UniProt ID #P31749

Storage: Supplied in 20 mM potassium phosphate (pH 7.0), 50 mM NaCl, 0.1 mM EDTA, 1 mg/ml BSA, 5% glycerol and 1% DMSO. Store at -20°C.

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.