

PAX9 (D9F1N) Rabbit mAb

Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Applications: W	Reactivity: H M R	Sensitivity: Endogenous	MW (kDa): 38	Source/Isotype: Rabbit IgG	UniProt ID: #P55771	Entrez-Gene Id: 5083
---------------------------	-----------------------------	-----------------------------------	------------------------	--------------------------------------	-------------------------------	--------------------------------

Product Usage Information**Application**

Western Blotting

Dilution

1:1000

Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

Specificity/Sensitivity

PAX9 (D9F1N) Rabbit mAb recognizes endogenous levels of total PAX9 protein.

Species predicted to react based on 100% sequence homology

Monkey

Source / Purification

Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Asp141 of human PAX9 protein.

Background

Paired box (PAX) proteins are a family of transcription factors that play important and diverse roles in animal development (1). Nine PAX proteins (PAX1-9) have been described in humans and other mammals. They are defined by the presence of an amino-terminal "paired" domain, consisting of two helix-turn-helix motifs, with DNA binding activity (2). PAX proteins are classified into four structurally distinct subgroups (I-IV) based on the absence or presence of a carboxy-terminal homeodomain and a central octapeptide region. Subgroup I (PAX1 and 9) contains the octapeptide but lacks the homeodomain; subgroup II (PAX2, 5, and 8) contains the octapeptide and a truncated homeodomain; subgroup III (PAX3 and 7) contains the octapeptide and a complete homeodomain; and subgroup IV (PAX4 and 6) contains a complete homeodomain but lacks the octapeptide region (2). PAX proteins play critically important roles in development by regulating transcriptional networks responsible for embryonic patterning and organogenesis (3); a subset of PAX proteins also maintain functional importance during postnatal development (4). Research studies have implicated genetic mutations that result in aberrant expression of PAX genes in a number of cancer subtypes (1-3), with members of subgroups II and III identified as potential mediators of tumor progression (2). PAX9 is expressed in pharyngeal arch mesenchyme (5,6), and is essential for embryonic development of the teeth and other pharyngeal arch derivatives (5,7). In mice, deletion of PAX9 results in an absence of structures derived from the pharyngeal pouches (e.g. thymus, parathyroid glands) (5), while in humans, PAX9 mutations are frequently associated with congenital tooth agenesis syndromes such as oligodontia and hypodontia (8,9). PAX9 appears to interact with other transcription factors (e.g. MSX1) to regulate the expression of BMP-family proteins (e.g. BMP4) that orchestrate pharyngeal arch development (10,11). Increased expression of PAX9, resulting from amplification at the PAX9 locus, has also been linked in research studies to an increased incidence of lung cancer (12).

Background References

- Lang, D. et al. (2007) *Biochem Pharmacol* 73, 1-14.
- Robson, E.J. et al. (2006) *Nat Rev Cancer* 6, 52-62.
- Wang, Q. et al. (2008) *J Cell Mol Med* 12, 2281-94.
- Blake, J.A. et al. (2008) *Dev Dyn* 237, 2791-803.
- Peters, H. et al. (1998) *Genes Dev* 12, 2735-47.
- Peters, H. et al. (1998) *Eur J Oral Sci* 106 Suppl 1, 38-43.
- Hetzer-Egger, C. et al. (2002) *Eur J Immunol* 32, 1175-81.
- Bergendal, B. et al. (2011) *Am J Med Genet A* 155A, 1616-22.
- Suda, N. et al. (2011) *J Dent Res* 90, 382-6.
- Nakatomi, M. et al. (2010) *Dev Biol* 340, 438-49.
- Wang, Y. et al. (2011) *J Dent Res* 90, 311-6.
- Kendall, J. et al. (2007) *Proc Natl Acad Sci U S A* 104, 16663-8.

Species Reactivity

Species reactivity is determined by testing in at least one approved application (e.g., western blot).

Western Blot Buffer

IMPORTANT: For western blots, incubate membrane with diluted primary antibody in 5% w/v nonfat dry milk, 1X TBS, 0.1% Tween® 20 at 4°C with gentle shaking, overnight.

Applications Key

W: Western Blotting

Cross-Reactivity Key

H: Human **M:** Mouse **R:** Rat

Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.