

#12931 Store at Room Temperature

Nuclease-free Water

50 ml



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 11/07/17

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Ultrapure Nuclease-free Water from Cell Signaling Technology is non-DEPC treated, endotoxin-tested, and contains no detectable endo- and exo-nucleases. Following stringent purification, it is 0.1 µm filtered into a polycarbonate bottle and autoclaved. It is suitable for molecular biology applications that require water devoid of nucleases.

Molecular Formula: H₂O

Molecular Weight: 18.02 g/mol

Storage: Store at room temperature. This product is stable for 24 months.

Directions for Use: Nuclease-free water is ready to use. Maintain sterility by performing aseptic technique when handling.

Please visit www.cellsignal.com for a complete listing of recommended companion products.