

Store at
RT

Water, LC-MS Grade (Burdick and Jackson™)

#27732

1L



Cell Signaling
TECHNOLOGY®

Support: +1-978-867-2388 (U.S.)
www.cellsignal.com/support

Orders: 877-616-2355 (U.S.)
orders@cellsignal.com

rev. 12/04/18

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Water, LC-MS Grade (Burdick and Jackson™) is highly purified and recommended for use in HPLC, UPLC, and mass spectrometry applications (LC-MS) for low background contaminants.

Storage: This product is stable for 18 months when stored at room temperature.

Burdick and Jackson is a trademark of Honeywell International, Inc.

Thank you for your recent purchase. If you would like to provide a review visit cellsignal.com/comments.

www.cellsignal.com

© 2018 Cell Signaling Technology, Inc.

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide **Species Cross-Reactivity:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.