

Store at
4°C

FoxP3/Transcription Factor Fixation/ Permeabilization Concentrate (4X)



Cell Signaling
TECHNOLOGY®

#44931

50 ml

Support: +1-978-867-2388 (U.S.)
www.cellsignal.com/support

Orders: 877-616-2355 (U.S.)
orders@cellsignal.com

New 02/19

For Research Use Only. Not For Use In Diagnostic Procedures.

Background: This product is a concentrated fixation and permeabilization reagent designed for use with components of the FoxP3/Transcription Factor Fixation/Permeabilization Kit #43481. This reagent immobilizes cellular proteins and enables antibody binding of transcription factors localized within the nucleus, and is compatible with fluorescent detection by flow cytometry.

Note: Precipitation may occur. The presence of precipitate does not affect the performance of the reagent.

Storage: Store at 4°C. Stable for at least 6 months when stored at 4°C.

Directions for Use: Dilute FoxP3/Transcription Factor Fixation/Permeabilization Concentrate (4X) #44931 to a 1X working concentration with FoxP3/Transcription Factor Fixation/Permeabilization Diluent (1X) #58766 prior to use. Combine 1 part Concentrate (4X) #44931 with 3 parts Diluent (1X) #58766. See the protocol for the FoxP3/Transcription Factor Fixation/Permeabilization Kit #43481 for further details.

For product specific protocols and a complete listing of recommended companion products please see the product web page at www.cellsignal.com.

Thank you for your recent purchase. If you would like to provide a review visit cellsignal.com/comments.

www.cellsignal.com

© 2019 Cell Signaling Technology, Inc.

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide **Species Cross-Reactivity:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.