



Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

Store at -20C
#4871

HSP40 (C64B4) Rabbit mAb

For Research Use Only. Not for Use in Diagnostic Procedures.

Applications: W, IHC-P	Reactivity: H M R Mk	Sensitivity: Endogenous	MW (kDa): 40	Source/Isotype: Rabbit	UniProt ID: #P25685	Entrez-Gene Id: 3337
----------------------------------	--------------------------------	-----------------------------------	------------------------	----------------------------------	-------------------------------	--------------------------------

Product Usage Information

Application

Western Blotting
Immunohistochemistry (Paraffin)

Dilution

1:1000
1:50

Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

Specificity/Sensitivity

HSP40 (C64B4) Rabbit mAb detects endogenous levels of total HSP40 protein. This antibody does not cross-react with other HSPs.

Source / Purification

Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to Glu223 of human HSP40/Hdj1.

Background

HSP40 and HSP40-like proteins represent a large family of chaperone proteins that are homologous to *E. coli* DnaJ protein (1). These proteins are classified into three subtypes based on their structures. The common feature of the family is a conserved J domain, which is usually located at the amino terminus of proteins and responsible for their association with HSP70 (1,2). Human HSP40, also known as Hdj1, belongs to subtype II that contain a unique Gly/Phe-rich region (2). HSP40 family proteins bind unfolded proteins, prevent their aggregation, and then deliver them to HSP70 (2,3). Another major function of HSP40 is to stimulate ATPase activity of HSP70, which causes conformational change of the unfolded proteins (4,5). The HSP40-HSP70-unfolded protein complex further binds to co-chaperones Hip, Hop and HSP90 or components of the protein degradation machinery such as CHIP and BAG-1, which either leads to protein folding or degradation, respectively (6).

Background References

- Cheetham, M.E. and Caplan, A.J. (1998) *Cell Stress Chaperones* 3, 28-36.
- Fan, C.Y. et al. (2003) *Cell Stress Chaperones* 8, 309-316.
- Langer, T. et al. (1992) *Nature* 356, 683-689.
- Liberek, K. et al. (1991) *Proc. Natl. Acad. Sci. USA* 88, 2874-2878.
- Cyr, D.M. et al. (1992) *J. Biol. Chem.* 267, 20927-20931.
- Höhfeld, J. et al. (2001) *EMBO Rep.* 2, 885-890.

Species Reactivity

Species reactivity is determined by testing in at least one approved application (e.g., western blot).

Western Blot Buffer

IMPORTANT: For western blots, incubate membrane with diluted primary antibody in 5% w/v BSA, 1X TBS, 0.1% Tween® 20 at 4°C with gentle shaking, overnight.

Applications Key

W: Western Blotting **IHC-P:** Immunohistochemistry (Paraffin)

Cross-Reactivity Key

H: Human **M:** Mouse **R:** Rat **Mk:** Monkey

Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

U.S. Patent No. 7,429,487, foreign equivalents, and child patents deriving therefrom.

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.