

Store at
-20°C

#50748

Astrocyte Markers Antibody Sampler Kit

Support: +1-978-867-2388 (U.S.)
www.cellsignal.com/supportOrders: 877-616-2355 (U.S.)
orders@cellsignal.com

New 01/21

For Research Use Only. Not For Use In Diagnostic Procedures.

Products Included	Product #	Quantity	Mol. Wt.	Isotype/Source
GFAP (E4L7M) XP® Rabbit mAb	80788	20 µl	50 kDa	Rabbit IgG
EAAT1 (D44E2) XP® Rabbit mAb	5684	20 µl	58 kDa	Rabbit IgG
Survivin (71G4B7) Rabbit mAb	2808	20 µl	16 kDa	Rabbit IgG
EAAT2 (E3P5K) Rabbit mAb	20848	20 µl	65 kDa	Rabbit IgG
ALDH1L1 (E7I2Q) Rabbit mAb	85828	20 µl	98 kDa	Rabbit IgG
AQP4 (D1F8E) XP® Rabbit mAb	59678	20 µl	28 kDa	Rabbit IgG
GAT1 (E7J1B) Rabbit mAb	37342	20 µl	65 kDa	Rabbit IgG
Anti-rabbit IgG, HRP-linked Antibody	7074	100 µl		Goat

See www.cellsignal.com for individual component applications, species cross-reactivity, dilutions, and additional application protocols.

Description: The Astrocyte Markers Antibody Sampler Kit provides an economical means of detecting astrocyte markers and control antibodies. The kit includes enough antibodies to perform two western blot experiments with each primary antibody.

Background: Astrocytes are a population of cells with distinctive morphological and functional characteristics that differ within specific areas of the brain. Postnatally, astrocyte progenitors migrate to reach their brain area and related properties. They have a regulatory role in brain functions that are implicated in neurogenesis and synaptogenesis, controlling blood-brain barrier permeability and maintaining extracellular homeostasis. Mature astrocytes also express some genes enriched in cell progenitors, suggesting they can retain proliferative potential (1). Astrocytes in the human brain are characterized by a heavy expression of the glial fibrillary acidic protein (GFAP), comprised of interlaminar that are located in layers I and II, protoplasmic in layers III and IV, varicose projections in layers V and VI, and fibrous astroglia in white matter (2,3).

Aquaporins (AQP) are integral membrane proteins that serve as channels in the transfer of water and small solutes across the membrane. AQP4 is present in the brain and it is enriched in astrocytes to regulate water homeostasis, preventing cerebral edema caused by solute imbalance (4,5). Excitatory amino acid transporters are members of sodium-dependent, high-affinity transporters that mediate the uptake of L-glutamate and D-aspartate (6). GABA transmitters are Na⁺/Cl⁻-dependent transporters that regulate neurotransmitter transport, including GAT1 (SLC6), GAT2, GAT3, and BGT1 (7). GAT1 expresses in the brain, preferentially in glial cells, but is also found in neurons, regulating the uptake and release of neurotransmitters in terminal clefts (8-10). EAAT2, also known as GLT-1, is primarily expressed in astrocytes accounting for up to 90% of the total glutamate transport in the brain. EAAT1 upregulates increased concentrations of glutamate in astrocytes and has a neuroprotective potential following ischemia since reactive astrocytes and activated microglia express EAAT1 but not EAAT2 (11-13).

Aldehyde dehydrogenase 1 family member L1 (ALDH1L1) is a member of the aldehyde dehydrogenase super family (14,15). ALDH1L1 has also been shown to be a useful astrocyte marker throughout the grey and white matter of the brain, labeling both the cell body and processes of astrocytes. ALDH1L1 does not label neurons or oligodendrocytes (16). It has been shown that expression of ALDH1L1 in astrocytes is stably expressed through normal tissue and during astrogliosis (17). Survivin is a 16 kDa anti-apoptotic protein highly expressed during fetal development and cancer malignancy. It binds and inhibits caspase-3, controlling the checkpoint in the G2/M-phase of the cell cycle by inhibiting apoptosis and promoting cell division (18). Survivin expression is associated with malignant phenotypes and prognosis of glioma (19).

Specificity/Sensitivity: Each antibody in the Astrocyte Markers Antibody Sampler Kit detects endogenous levels of its target protein. ALDH1L1 (E7I2Q) Rabbit mAb recognizes endogenous levels of total ALDH1L1 protein. This antibody detects a non-specific band of unknown origin in rodent tissue extracts at 47 kDa.

Source/Purification: Monoclonal antibodies are produced by immunizing rabbits with synthetic peptides corresponding to residues near the amino terminus of human GFAP, EAAT2, and EAAT1 proteins. Monoclonal antibodies are produced by immunizing rabbits with synthetic peptides corresponding to residues near the carboxy terminus of human GAT1 protein and specific to the carboxy terminus of human AQP4 protein. Monoclonal antibodies are produced by immunizing rabbits with synthetic peptides corresponding to residues surrounding Val431 of human ALDH1L1 protein and Cys60 of human Survivin protein.

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibodies.

Please visit www.cellsignal.com for validation data and a complete listing of recommended companion products.

Background References:

- (1) Siracusa, R. et al. (2019) *Front Pharmacol* 10, 1114.
- (2) Vasile, F. et al. (2017) *Brain Struct Funct* 222, 2017-2029.
- (3) Eng, L.F. et al. (2000) *Neurochem Res* 25, 1439-51.
- (4) Takata, K. et al. (2004) *Prog Histochem Cytochem* 39, 1-83.
- (5) Kobayashi, H. et al. (2004) *J Pharmacol Sci* 96, 264-70.
- (6) Amara, S.G. and Fontana, A.C. (2002) *Neurochem Int* 41, 313-8.
- (7) Kristensen, A.S. et al. (2011) *Pharmacol Rev* 63, 585-640.
- (8) Borden, L.A. (1996) *Neurochem Int* 29, 335-56.
- (9) Moldavan, M. et al. (2017) *J Neurophysiol* 118, 3092-3106.
- (10) Lorenz-Guertin, J.M. and Jacob, T.C. (2018) *Dev Neurobiol* 78, 238-270.
- (11) Hediger, M.A. (1999) *Am J Physiol* 277, F487-92.
- (12) Gegelashvili, G. et al. (1996) *Neuroreport* 8, 261-5.
- (13) Beschorner, R. et al. (2007) *Histopathology* 50, 897-910.
- (14) Krupenko, S.A. (2009) *Chem Biol Interact* 178, 84-93.
- (15) Krupenko, S.A. and Oleinik, N.V. (2002) *Cell Growth Differ* 13, 227-36.
- (16) Cahoy, J.D. et al. (2008) *J Neurosci* 28, 264-78.
- (17) Michalovicz, L.T. et al. (2019) *J Neurochem* 150, 420-440.
- (18) Reed, J.C. and Reed, S.I. (1999) *Nat Cell Biol* 1, E199-200.
- (19) Tong, X. et al. (2019) *Oncol Lett* 18, 359-367.

Thank you for your recent purchase. If you would like to provide a review visit www.cellsignal.com/comments.

www.cellsignal.com

© 2021 Cell Signaling Technology, Inc.

XP and Cell Signaling Technology are trademarks of Cell Signaling Technology, Inc.

Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide Species Cross-Reactivity: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.