

Store at
4°C

#58766

FoxP3/Transcription Factor Fixation/ Permeabilization Diluent (1X)

160 ml



Cell Signaling
TECHNOLOGY®

Support: +1-978-867-2388 (U.S.)
www.cellsignal.com/support

Orders: 877-616-2355 (U.S.)
orders@cellsignal.com

New 01/19

For Research Use Only. Not For Use In Diagnostic Procedures.

Background: This product is a diluent designed for use with components of the FoxP3/Transcription Factor Fixation/Permeabilization Kit #43481. Proper dilution of the FoxP3/Transcription Factor Fixation/Permeabilization Concentrate (4X) #44931 results in the reagent immobilizing cellular proteins and enabling antibody binding of transcription factors localized within the nucleus. This product is compatible with fluorescent detection by flow cytometry.

Note: Precipitation may occur. The presence of precipitate does not affect the performance of the reagent.

Storage: Store at 4°C. Stable for 6 months when stored at 4°C.

Directions for Use: Use FoxP3/Transcription Factor Fixation/Permeabilization Diluent (1X) #58766 to dilute FoxP3/Transcription Factor Fixation/Permeabilization Concentrate (4X) #44931 to a 1X working concentration prior to use. Combine 3 parts Diluent (1X) #58766 with 1 part Concentrate (4X) #44931. See the protocol for the FoxP3/Transcription Factor Fixation/Permeabilization Kit #43481 for further details.

For product specific protocols and a complete listing of recommended companion products please see the product web page at www.cellsignal.com.

Thank you for your recent purchase. If you would like to provide a review visit www.cellsignal.com/comments.

www.cellsignal.com

© 2019 Cell Signaling Technology, Inc.

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide **Species Cross-Reactivity:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.