

#6532 Store at -20°C

# SignalSilence® MEK2 siRNA II



✓ 10 µM in 300 µl  
(100 Transfections)

**Orders** ■ 877-616-CELL (2355)  
orders@cellsignal.com  
**Support** ■ 877-678-TECH (8324)  
info@cellsignal.com  
**Web** ■ www.cellsignal.com

rev. 02/11/16

**For Research Use Only. Not For Use In Diagnostic Procedures.**


### Species Cross-Reactivity: H, M, R

**Description:** SignalSilence® MEK2 siRNA from Cell Signaling Technology allows the researcher to specifically inhibit MEK2 expression using RNA interference, a method in which gene expression can be selectively silenced through the delivery of double stranded RNA molecules into the cell. All SignalSilence® siRNA products are rigorously tested in-house and have been shown to reduce target protein expression by western analysis.

**Background:** MEK1 and MEK2, also called MAPK or Erk kinases, are dual-specificity protein kinases that function in a mitogen activated protein kinase cascade controlling cell growth and differentiation (1-3). Activation of MEK1 and MEK2 occurs through phosphorylation of two serine residues at positions 217 and 221 (in the activation loop of subdomain VIII) by Raf-like molecules. MEK1/2 is activated by a wide variety of growth factors and cytokines and also by membrane depolarization and calcium influx (1-4). Constitutively active forms of MEK1/2 are sufficient for the transformation of NIH/3T3 cells or the differentiation of PC12 cells (4). MEK activates p44 and p42 MAP kinase by phosphorylating both threonine and tyrosine residues at sites located within the activation loop of kinase subdomain VIII.

**Directions for Use:** CST recommends transfection with 100 nM MEK2 siRNA 48 to 72 hours prior to cell lysis. For transfection procedure, follow protocol provided by the transfection reagent manufacturer. Please feel free to contact CST with any questions on use.

**Quality Control:** Oligonucleotide synthesis is monitored base by base through trityl analysis to ensure appropriate coupling efficiency. The oligo is subsequently purified by affinity-solid phase extraction. The annealed RNA duplex is further analyzed by mass spectrometry to verify the exact composition of the duplex. Each lot is compared to the previous lot by mass spectrometry to ensure maximum lot-to-lot consistency.



Western blot analysis of extracts from HeLa cells, transfected with 100 nM SignalSilence® Control siRNA (Fluorescein Conjugate) #6201 (-) or SignalSilence® MEK2 siRNA II (+), using MEK2 (13E3) Rabbit mAb #9147 and α-Tubulin (11H10) Rabbit mAb #2125. MEK2 (13E3) Rabbit mAb confirms silencing of MEK2 expression, while the α-Tubulin (11H10) Rabbit mAb is used to control for loading and specificity of MEK2 siRNA.

Entrez-Gene ID #5605  
Swiss-Prot Acc. #P36507

**Storage:** MEK2 siRNA II is supplied in RNase-free water. Aliquot and store at -20°C.

**Please visit [www.cellsignal.com](http://www.cellsignal.com) for a complete listing of recommended companion products.**

### Background References:

- (1) Crews, C.M. et al. (1992) *Science* 258, 478-480.
- (2) Alessi, D.R. et al. (1994) *EMBO J.* 13, 1610-1619.
- (3) Rosen, L.B. et al. (1994) *Neuron* 12, 1207-1221.
- (4) Cowley, S. et al. (1994) *Cell* 77, 841-852.

© 2009 Cell Signaling Technology, Inc.

**Applications Key:** W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide  
**Species Cross-Reactivity Key:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine  
Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.