

Ubiquityl-Histone H2A.Z (Lys120) (E6G7X) Rabbit mAb



Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Applications:	Reactivity:	Sensitivity:	MW (kDa):	Source/Isotype:	UniProt ID:	Entrez-Gene Id:
W, ChIP, ChIP-seq	H M R Mk	Endogenous	23	Rabbit IgG	#P0C0S5	3015

Product Usage Information

For optimal ChIP and ChIP-seq results, use 10 µl of antibody and 10 µg of chromatin (approximately 4×10^6 cells) per IP. This antibody has been validated using SimpleChIP® Enzymatic Chromatin IP Kits.

Application

Western Blotting
Chromatin IP
Chromatin IP-seq

Dilution

1:1000
1:50
1:50

Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. *Do not aliquot the antibody.*

Specificity/Sensitivity

Ubiquityl-Histone H2A.Z (Lys120) (E6G7X) Rabbit mAb recognizes endogenous levels of histone H2A.Z protein only when ubiquitylated at Lys120. This antibody does not detect histone H2A.Z that is ubiquitylated on Lys121, nor does it detect histone H2A.Z that is dually ubiquitylated on Lys120 and Lys121. This antibody does not cross-react with other ubiquitylated proteins or free ubiquitin.

Source / Purification

Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues surrounding ubiquitylated Lys120 of human histone H2A.Z protein.

Background

Modulation of chromatin structure plays a critical role in the regulation of transcription in eukaryotes. The nucleosome, made up of four core histone proteins (H2A, H2B, H3 and H4), is the primary building block of chromatin. In addition to the growing number of post-translational histone modifications regulating chromatin structure, cells can also exchange canonical histones with variant histones that can directly or indirectly modulate chromatin structure (1). There are five major variants of histone H2A: canonical H2A (most abundant), H2A.X, MacroH2A, H2ABbd and H2A.Z (2). Histone H2A.Z, the most conserved variant across species, functions as both a positive and negative regulator of transcription and is important for chromosome stability (2). Several homologous protein complexes, such as SWR-C (*S. cerevisiae*), TIP60 (*D. melanogaster*) and SRCAP (mammals), have been shown to catalyze the ATP-dependent exchange of H2A.Z for H2A in the nucleosome (3,4,5). This exchange of histone H2A variants changes histone-histone interactions in the nucleosome core and alters an acidic patch on the surface of the nucleosome, resulting in changes in nucleosome stability and binding of non-histone proteins such as HP1α (6,7).

H2A.Z is a histone H2A variant protein that is critical for proper regulation of gene expression. H2A.Z is localized throughout the genome, but appears to be most concentrated at the promoters and enhancers of active genes (8). Acetylation of histone H2A.Z at promoters and enhancers confers nucleosome destabilization and open chromatin confirmation, facilitating transcriptional activation (9-11). While the bulk of histone H2A.Z appears to be excluded from constitutive heterochromatin, histone H2A.Z is found in various forms of facultative heterochromatin, including the inactive X chromosome and transcriptionally poised bivalent gene promoters (8,12). In these heterochromatic regions of the genome, H2A.Z is mono-ubiquitylated on Lys120 and Lys121 by the Ring1B ubiquitin ligase found in the Polycomb Repressor Complex 1 (PRC1) (8,12). Mono-methylation of H2A.Z on Lys120 and Lys121 facilitates repression of gene expression by inhibiting the binding of the activating BRD4 protein and facilitating the recruitment of the Polycomb Repressor Complex 2 (PRC2), the latter of which methylates histone H3 on Lys27 and facilitates transcriptional repression (8).

Background References

- Jin, J. et al. (2005) *Trends Biochem Sci* 30, 680-7.
- Raisner, R.M. and Madhani, H.D. (2006) *Curr Opin Genet Dev* 16, 119-24.
- Mizuguchi, G. et al. (2004) *Science* 303, 343-8.
- Kusch, T. et al. (2004) *Science* 306, 2084-7.
- Ruhl, D.D. et al. (2006) *Biochemistry* 45, 5671-7.
- Suto, R.K. et al. (2000) *Nat Struct Biol* 7, 1121-4.
- Fan, J.Y. et al. (2004) *Mol Cell* 16, 655-61.
- Surface, L.E. et al. (2016) *Cell Rep* 14, 1142-55.
- Millar, C.B. et al. (2006) *Genes Dev* 20, 711-22.

10. Ishibashi, T. et al. (2009) *Biochemistry* 48, 5007-17.
11. Valdés-Mora, F. et al. (2012) *Genome Res* 22, 307-21.
12. Sarcinella, E. et al. (2007) *Mol Cell Biol* 27, 6457-68.
-

Species Reactivity	Species reactivity is determined by testing in at least one approved application (e.g., western blot).
Western Blot Buffer	IMPORTANT: For western blots, incubate membrane with diluted primary antibody in 5% w/v BSA, 1X TBS, 0.1% Tween® 20 at 4°C with gentle shaking, overnight.
Applications Key	W: Western Blotting ChIP: Chromatin IP ChIP-seq: Chromatin IP-seq
Cross-Reactivity Key	H: Human M: Mouse R: Rat Mk: Monkey
Trademarks and Patents	Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc. SimpleChIP is a registered trademark of Cell Signaling Technology, Inc. All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.
Limited Uses	Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect. Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.