

# $\beta$ -Actin (13E5) Rabbit mAb (Alexa Fluor® 594 Conjugate)



**Orders** ■ 877-616-CELL (2355)  
orders@cellsignal.com

**Support** ■ 877-678-TECH (8324)  
info@cellsignal.com

**Web** ■ www.cellsignal.com

rev. 04/14/16

**For Research Use Only. Not For Use In Diagnostic Procedures.**

Applications	Species Cross-Reactivity*	Isotype
IF-IC Endogenous	H, M, R, Mk, B, Pg, (C, Dg, Hr)	Rabbit IgG


**Description:** This Cell Signaling Technology antibody is conjugated to Alexa Fluor® 594 fluorescent dye and tested in-house for immunofluorescent analysis in monkey cells. This antibody is expected to exhibit the same species cross-reactivity as the unconjugated  $\beta$ -Actin (13E5) Rabbit mAb #4970.

**Background:** Actin, a ubiquitous eukaryotic protein, is the major component of the cytoskeleton. At least six isoforms are known in mammals. Nonmuscle  $\beta$ - and  $\gamma$ -actin, also known as cytoplasmic actin, are predominantly expressed in nonmuscle cells, controlling cell structure and motility (1).  $\alpha$ -cardiac and  $\alpha$ -skeletal actin are expressed in striated cardiac and skeletal muscles, respectively; two smooth muscle actins,  $\alpha$ - and  $\gamma$ -actin, are found primarily in vascular smooth muscle and enteric smooth muscle, respectively. These actin isoforms regulate the contractile potential of muscle cells (1). Actin exists mainly as a fibrous polymer, F-actin. In response to cytoskeletal reorganizing signals during processes such as cytokinesis, endocytosis, or stress, cofilin promotes fragmentation and depolymerization of F-actin, resulting in an increase in the monomeric globular form, G-actin (2). The Arp2/3 complex stabilizes F-actin fragments and promotes formation of new actin filaments (2). It has been reported that actin is hyperphosphorylated in primary breast tumors (3). Cleavage of actin under apoptotic conditions has been observed *in vitro* and in cardiac and skeletal muscle (4-6). Actin cleavage by caspase-3 may accelerate ubiquitin/proteasome-dependent muscle proteolysis (6).

**Specificity/Sensitivity:**  $\beta$ -Actin (13E5) Rabbit mAb (Alexa Fluor® 594 Conjugate) detects endogenous levels of total  $\beta$ -actin protein. This antibody may cross-react with the  $\gamma$ -actin (cytoplasmic isoform). It does not cross-react with  $\alpha$ -skeletal,  $\alpha$ -cardiac,  $\alpha$ -vascular smooth, or  $\gamma$ -enteric smooth muscle isoforms.

**Source/Purification:** Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues near the amino terminus of human  $\beta$ -actin protein.

## COS-7



Confocal immunofluorescent analysis of COS-7 cells using  $\beta$ -Actin (13E5) Rabbit mAb (Alexa Fluor® 594 Conjugate) (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye).

Entrez-Gene ID #60  
UniProt ID #P60709

**Storage:** Supplied in PBS (pH 7.2), less than 0.1% sodium azide and 2 mg/ml BSA. Store at 4°C. *Protect from light. Do not freeze.*

**\*Species cross-reactivity other than monkey is determined by western blot using the unconjugated antibody.**

### Recommended Antibody Dilutions:

Immunofluorescence (IF-IC) 1:50  
IF Protocol: Methanol Permeabilization required

**For application specific protocols please see the web page for this product at [www.cellsignal.com](http://www.cellsignal.com).**

**Please visit [www.cellsignal.com](http://www.cellsignal.com) for a complete listing of recommended companion products.**

### Background References:

- (1) Herman, I.M. (1993) *Curr. Opin. Cell Biol.* 5, 48-55.
- (2) Condeelis, J. (2001) *Trends Cell Biol.* 11, 288-293.
- (3) Lim, Y.P. et al. (2004) *Clin. Cancer Res.* 10, 3980-3987.
- (4) Kayalar, C. et al. (1996) *Proc. Natl. Acad. Sci. USA.* 93, 2234-2238.
- (5) Communal, C. et al. (2002) *Proc. Natl. Acad. Sci. USA.* 99, 6252-6256.
- (6) Du, J. et al. (2004) *J. Clin. Invest.* 113, 115-123.

U.S. Patent No. 5,675,063

Alexa Fluor® is a registered trademark of Molecular Probes, Inc.

DRAQ5® is a registered trademark of Biostatus Limited.