

#9871
Store at Room
Temperature

20 mM Citrate pH 3.0 (Sterile)

- Small 10 ml
 Large 50 ml



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com

rev. 01/05/21

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: 20 mM Citrate pH 3.0 (Sterile) is specifically formulated to be used as a buffer for reconstituting sterile cytokines or antibodies.

Endotoxin: Less than 0.5 EU/ml citrate buffer.

Storage: Store at room temperature. Maintain sterility.