

Store at Room
Temperature

#9999

Nonfat Dry Milk

- 250 grams
- Sample Size



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 11/07/17

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Premium quality nonfat dry milk for use as a blocking agent during Western blotting to block non-specific antibody binding and reduce background signal.

Storage: Store at room temperature.

Companion Products:

Tris Buffered Saline with Tween 20 (TBST - 10X) #9997

BSA #9998

Prestained Protein Marker, Broad Range (Premixed Format) #7720

Biotinylated Protein Ladder Detection Pack #7727

20X LumiGLO® Reagent and 20X Peroxide #7003

Phototope®-HRP Western Blot Detection System, Anti-rabbit IgG, HRP-linked Antibody #7071

Phototope®-HRP Western Blot Detection System, Anti-mouse IgG, HRP-linked Antibody #7072

Anti-rabbit IgG, HRP-linked Antibody #7074

Anti-biotin, HRP-linked Antibody #7075

Anti-mouse IgG, HRP-linked Antibody #7076

Anti-rat IgG, HRP-linked Antibody #7077