



**Phospho-Akt (Ser473) (D9E) XP® Rabbit mAb #4060** (upper) and **Akt (pan) (C67E7) Rabbit mAb #4691** (lower).  
For full validation details see [www.cellsignal.com](http://www.cellsignal.com)

The Highest Quality

# Akt Antibodies

are from Cell Signaling Technology

**Over 45 thoroughly validated Akt Antibodies available**

Cell Signaling Technology offers the widest range of phospho and total protein Akt antibodies. These reagents are the most thoroughly validated, widely used, and highly cited Akt antibodies available.

- :: CST™ phospho and total protein Akt antibodies offer **unsurpassed specificity, sensitivity, and performance.**
- :: **Extensive in-house validation** in a wide range of applications means that optimization is not left up to you, the user.
- :: **Technical support is provided by scientists** who produce the products and know them best.

Unparalleled product quality, validation, and technical support

for quality products you can trust...

[www.cellsignal.com](http://www.cellsignal.com)

