



Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) (D13.14.4E) XP® Rabbit mAb #4370 (top) and p44/42 MAPK (Erk1/2) (137F5) Rabbit mAb #4695 (bottom).

For full validation details see www.cellsignal.com

The Highest Quality MAPK Antibodies

are from Cell Signaling Technology

Unparalleled product quality, validation, and technical support

Over 30 thoroughly validated MAPK Antibodies available

Cell Signaling Technology (CST™) offers the widest range of phospho and total protein MAPK antibodies. These reagents are the most thoroughly validated, widely used, and highly cited MAPK antibodies available.

- ❑ CST™ phospho and total protein MAPK antibodies offer **unsurpassed specificity, sensitivity, and performance.**
- ❑ **Extensive in-house validation** in a wide range of applications means that optimization is not left up to you, the user.
- ❑ **Technical support is provided by scientists** who produce the products and know them best.



WINNER! Life Science Industry Awards®
 ❑ Best Antibodies ❑ Best Breakthrough Products for Cancer Research

for quality products you can trust...
www.cellsignal.com

