

The Highest Quality Cytokines & Growth Factors


...from Cell Signaling Technology

The world's highest quality antibody provider has now extended its expertise to cytokine production.

- ⚡ Produced and bioassayed in house with the highest purity and bioactivity.
- ⚡ Comparison of multiple lots, stringent product specifications, and rigorous quality control ensures maximum lot-to-lot consistency.
- ⚡ Endotoxin levels are routinely tested and are less than 0.01 ng/μg cytokine.
- ⚡ Many products are produced in mammalian cells to maximize natural conformation and glycosylation.
- ⚡ Multi-milligram quantities available.


Human Interleukin-17A (hIL-17A) #8928

Purity




6 μg reduced (+) and non-reduced (-) recombinant hIL-17A.

Bioactivity



The production of IL-6 by human foreskin fibroblasts cultured with increasing concentrations of hIL-17A for 48 hours was assessed.

Downstream Signaling



Human foreskin fibroblasts treated with hIL-17A for 15 minutes were examined using phospho- and total p38 MAPK antibodies.

for quality products you can trust...

www.cellsignal.com

China www.cst-c.com.cn | Europe www.cellsignal.eu | Japan www.cstj.co.jp



Cell Signaling
TECHNOLOGY®

Tel: 978-867-2300 | E-mail: info@cellsignal.com